

AQUATEC®

Ocean Ergo Family

Clinical Introduction





Ocean Ergo
Family designed
to make personal
care activies
easier for both the
user and carer

Contents

Introduction	04
Meet the family	04
Seating	06
Creating stability	08
What is immersion?	10
Off-loading	12
Envelopment	12
Positioning for body systems	13
Positioning for waste elimination	14
The Puborectalis Muscle	15
The importance of the seat angle	15
Positioning and the Autonomic Nervous System	16
The nervous system in waste elimination	18
Ocean Accessories	20
Positioning accessories	21
Safety accessories	23
Seat and backrest accessories	24
References	27





Introduction to the **Ocean Ergo Family**

Personal care is an important part of a person's life and wellbeing. When personal care activities become difficult, a healthcare professional may recommend using a shower chair commode as part of a care plan.

The Aquatec Ocean Ergo Family of shower chair commodes have been designed with user comfort, independence and dignity in mind. Based on the well-established Aquatec Ocean range, the Ocean Ergo family comes packed with a range of new features, designed to make personal care activities easier for both the user and carer.

The purpose of this booklet is to give an overview of the importance of seating and accessory selection in relation to the Ocean Ergo Family of shower chair commodes.

For ease of reference, this guide has been split into **two** sections:

1 Seating 2 Accessories

The Aquatec Ocean Ergo Family of shower chair commodes

have been designed
with user comfort,
independence and
dignity in mind.

KEY FEATURES AT A GLANCE



Ergonomic seat

Insight-led design encourages a more upright seating position for enhanced independence and comfort.

Ergonomic tilt

Intuitive, ergonomic tilt with weight shifting abilities

Highly customisable

Can be tailored to each individual by making some simple adjustments or by adding some easy to fit, off the shelf accessories.

Meet the family

From a standard transit model, through to tilt and recline, meet the flexible range of shower chairs that offer personal care solutions for a wide range of physical conditions.

Ocean 24"
Ergo XL ▼

Ocean
Ergo XL ▼

Ocean VIP
Ergo ▼

Ocean
Ergo ▼

Ocean Dual
VIP Ergo ▼

Ocean 24"
Ergo ▼





Seating principles...


Don't change, no matter what seat surface you sit on!

2

- ▶ It's true, seating principles don't change, no matter the seat surface you are sitting on. Therefore, given how quickly pressure ulcer development can happen, it's important that every therapeutic surface offers a level of skin protection. The New Ocean ergonomic seat plate has been designed with this in mind. When seating a user on the new Ocean Ergo shower chair commodes, the aims of seating remain the same:

- 1 ▶ Normalise tone or decrease abnormal influence on the body.
- 2 ▶ Maintain skeletal alignment.
- 3 ▶ Prevent, accommodate or correct skeletal deformity.
- 4 ▶ Provide stable base of support to promote function.
- 5 ▶ Promote increased tolerance of desired position.
- 6 ▶ Promote comfort and relaxation.
- 7 ▶ Facilitate normal movement patterns or control abnormal movement patterns.
- 8 ▶ Manage pressure or prevent the development of pressure ulcers.
- 9 ▶ Decrease fatigue.
- 10 ▶ Enhance autonomic nervous system function (e.g. cardiac, digestive and respiratory function).
- 11 ▶ Facilitate maximum function with minimum pathology.

(Jones and Gray 2005)



Key to seating success

An interdependency between:

- ▶ Seating a user successfully is based on an interdependency between stability, comfort and function. Stability of the pelvis is key to ensuring overall comfort for the user, and only when the user is comfortable will they be able to relax and sit in a position which supports optimal functional performance in activities of daily living (ADL).



3

Creating stability

Stability is created via the pelvis, for this you need...

- Pelvic stability is dependent upon three things: immersion, off-loading and envelopment.

Immersion, off-loading and envelopment are generally thought of in relation to a seat cushion, so how are they created on a shower chair commode?

Immersion ▾



Off-loading ▾



Envelopment ▾





Bony landmarks you need to know

We need to also consider the difference:

In height between the IT's and the hip joints as this is important when it comes to ensuring immersion.

- ▶ To understand immersion, off-loading and envelopment it's important to be familiar with the anatomy of the pelvis, specifically the shape and location of the Ischial Tuberossities (IT's) and the hip joints. Important to note is the height difference between these two bony prominences, as this plays an important role when it comes to immersion and off-loading.

- Hip joints
- The Ischial Tuberossities (IT's)

What is Immersion?

Immersion is the depth to which a body penetrates (sinks into) a seat from an uppermost plane, for example, a seat cushion.

4


- ▶ The body can only do this if the material upon which it is seated permits it.
- ▶ The IT's are the first bony prominence to come into contact with any seat surface.
- ▶ Humans are not designed for sitting, our IT's are small and fairly pointed, which quickly leads to the build up of pressure and discomfort.
- ▶ As a result, we constantly change our position to relieve pressure and make ourselves more comfortable.

- Hip joints
- The Ischial Tuberossities (IT's)



Q. How is immersion created on a firmer surface such as the Ocean Ergonomic Seat Plate?

A. Via the seat aperture



- ▶ The aperture needs to be wide enough to allow the IT's to fit in it.
- ▶ The distance between the IT's on an adult human pelvis is roughly 12-14cm, and less in children.

Unlike on a seat cushion where the IT's can immerse into the foam or other material, on the Ocean Ergonomic Seat Plate, immersion of the IT's is created by allowing them to fit within the aperture of the seat plate.


Humans are all different shapes and sizes, however, the average distance between the IT's on an adult human pelvis, irrespective of how much adipose tissue someone has, is between 12-14cm, and less in children.

5

Off-loading

Off-loading is the clinical practice of reducing or removing pressure from one area of the body to another, in an effort to reduce risk of injury.

When the IT's immerse into the seat aperture, the hip joints come into contact with the seat surface, off-loading weight away from the IT's to these larger bony weight bearing surfaces. Without the seat aperture this would not happen and the user would be left balancing on their IT's, creating peak pressure build up under these bony prominences.



Weight transferred away from the IT's to the hip joints

6

Envelopment

Envelopment is the ability of a seat surface to conform, fit or mould around the irregular shape of the body.

The reality is that creating envelopment on a firmer surface is difficult. However, by incorporating a precontoured shape into the seat plate design based on the shape of an average adult pelvis, envelopment is improved as the contours of the seat plate respect the anatomical shape of the hips, trochanters and femurs, creating improved contact (envelopment) between user and seat.

- ▶ Envelopment is essential for pressure redistribution as it increases surface contact area between the user and the seat surface.
- ▶ Precontouring based on moulds taken of average adult pelvises respects the anatomical dimensions of the trochanters and femurs helping maximise surface contact area.
- ▶ This improves pressure redistribution, comfort, positioning and overall stability whilst seated.



$$\text{Pressure} = \frac{\text{FORCE}}{\text{AREA}}$$




7

Positioning is important for physiological body systems

The pelvis and the spine are connected therefore when we stabilise the pelvis, this encourages the spine to adopt its natural curvatures or s-shape, which extends the trunk, head and neck. This is essential, not only for optimal positioning, but also for physiological body system function.

- ▶ There are 11 physiological systems within the body, these include the digestive and the urinary system.
- ▶ Positioning is important when it comes to the overall function of these body systems.
- ▶ Slouching compresses the body's internal organs, including the digestive tract and the bladder, this can lead to problems with constipation, or urinary retention.
- ▶ When the pelvis is 'stable' the upper body (trunk), head and neck can extend. This upright sitting posture helps facilitate digestive and urinary system function, essential for overall health of the user.
- ▶ Ensuring an upright seated posture is key to the function of these and the other body systems, as a slouched posture will compress the body's internal organs which will affect their function.




8

Positioning for waste elimination

Positioning is also important when it comes to elimination of waste products from the body. Waste is stored in the colon, it then moves to the rectum and then the anus, where it is eliminated from the body.

Waste should be free to move through the colon and pass to the rectum and eventually the anus, where it is eliminated from the body. A slouched posture may compress the colon causing interruption to the function of this organ.



9


The Puborectalis Muscle

Positioning is important when it comes to the muscle which controls bowel continence, the puborectalis muscle.

Historically, humans adopted a squat position in order to eliminate their bowels. Over the centuries, however, more civilised ways were developed with the introduction of the toilet. The position adopted on a toilet, however, became more similar to that of sitting on a chair and the squat position became lost.

When examining the impact positioning had on the shape of the rectum and the puborectalis muscle, it was clear that the squat position was most effective in helping relax the puborectalis muscle, allowing the rectum to straighten and bowel elimination to become easier.

- ▶ When standing, the muscle tightens around the rectum, maintaining continence.
- ▶ When sitting it partially relaxes.
- ▶ When squatting, it fully relaxes, allowing the rectum to straighten, making elimination of waste easier.



10

The importance of the seat angle

In order to try and achieve more of a squat position, for more effective bowel elimination, the new ergonomic seat plate has been designed with a 5 degree rearward angle, which encourages the knees slightly closer to the chest helping to relax the puborectalis muscle.

5° rearward angle






Positioning and the Autonomic Nervous System (ANS)

11

The ANS which consists of the **Sympathetic Nervous System (SNS)** and the **Parasympathetic Nervous System (PNS)** is responsible for controlling the rectum and the bladder.

Going to the toilet is a private activity which we normally do in the comfort of our own bathrooms or WC's in our own homes, away from other people.

Consider, however, if you are seated on a shower chair commode and you are not positioned well. You may feel as though you are slipping down the chair or falling forward. If this is the case you're not likely to be relaxed and comfortable, which will prevent you from being able to eliminate your bowel or bladder, this is due to the functioning of the sympathetic and parasympathetic nervous systems.



12

The role of the body's nervous system in waste elimination

The PNS controls the body at rest, it is responsible for the body's "rest and digest" functions including contraction of the bladder and relaxation of the rectum, essential for waste elimination.

- ▶ The SNS on the other hand, controls the body's response to a perceived threat and is responsible for "fight or flight" reactions. When this system is activated, like in the case of feeling unstable or insecure whilst seated, the bladder relaxes and the rectum contracts, preventing waste elimination.
- ▶ The parasympathetic nervous system (PNS) controls homeostasis and the body at rest and is responsible for the body's "rest and digest" function. The sympathetic nervous system (SNS) controls the body's responses to a perceived threat and is responsible for the "fight or flight" response.



Homeostasis
**is a dynamic balance between
the autonomic branches**

- We can eliminate our bowel and bladder when the parasympathetic nervous system (our rest and digest system) is activated, which contracts the bladder and relaxes the rectum allowing us to urinate or pass faeces. If we are nervous or afraid, our sympathetic nervous system is triggered which relaxes the bladder and contracts the rectum, preventing us from passing urine or faeces.



► **Sympathetic**

● **Fight or flight:**

Sympathetic activity dominates

13

Ocean Accessories

Creating optimal positioning,
stability, comfort and function...

...but above all else, safety!



a

Positioning accessories **Lateral supports**

Lateral supports are available with all Ocean Ergo models.

- ▶ Delivered individually, however, two laterals should always be ordered and fitted to ensure optimal positioning and support for the user.
- ▶ Three points of height adjustment situated on the backrest cane.
- ▶ Additional height and depth adjustment on the lateral itself.
- ▶ The pad is broad and padded to ensure comfort and pressure redistribution when in use.

**b**

Positioning accessories **Headrest**

Comes as standard with all models (unless specified otherwise).

- ▶ Height and depth adjustable headrest.
- ▶ Height and depth are easily adjusted via the blue handscrews.
- ▶ Offers mild positioning and support.



c

Positioning accessories

Ergonomic headrest holder and head/neckrest

The ergonomic headrest holder comes as standard on all tilt in space models.

- ▶ The head/neckrest with ergonomic design, supports around the sub occipital head region.
- ▶ The ergonomic headrest holder offers various adjustment possibilities and is easily adjusted via the handscrews.
- ▶ Choice of headrest or head/neckrest with ergonomic head rest holder.

**d**

Positioning accessories

Calf and amputee support

The calf and amputee support offers both height and angle adjustment, with a contoured calf pad for added positioning and comfort.

- ▶ The calf and amputee support is also available with a foot support for users who may need their lower limbs to be elevated for medical reasons or comfort.
- ▶ All adjustments can be easily done by hand by loosening and tightening the handscrews.



e

Safety accessories

Pelvic strap and chest strap

The single pull pelvic strap offers additional security for the user, but can also be used for positioning purposes should the user have any postural tendencies.

- ▶ The pelvic strap is made from soft fabric and is broad to ensure comfort and reduce the risk of pressure build up.
- ▶ Positioning the belt between 45-60 degrees will help prevent sliding due to a posteriorly tilted pelvis, whilst attachment higher on the backrest canes will help to control the pelvis if anterior pelvic tilt is a problem.
- ▶ The strap has three points of adjustment and can be shortened or lengthened when the user is seated.
- ▶ The pelvic strap can also act as a chest strap in situations where anterior support is required. When using the chest strap, it should always be done in combination with the pelvic strap. It is essential that the pelvis is stable and well positioned before fitting any type of anterior support, to prevent risk of injury to the user.

**f**

Safety accessories

Anti-tip levers

Anti-tip levers can be fitted both at the rear and front of the Ocean Ergo models, for added stability and safety during transfer or when the chair is in tilt.

- ▶ In addition the anti-tip levers provide extra stability when used in showering situations where the floor surface is cambered for water drainage.



g**Safety accessories****Armrest locking device**

The armrest locking device is a safety feature which locks the armrest into a horizontal or vertical position.

- ▶ Especially useful on the Tilting Ocean models to prevent the armrests being lifted upwards unintentionally when the chair is tilted.
- ▶ Also useful if the user experiences uncontrolled spasms or movements.

**h****Seat and backrest accessories****Soft seat range**

Accessories such as the soft seat, which easily fits over the existing ergonomic seat plate, gives users with comfort, sensory or skin protection issues added relief and may help increase their sitting tolerance.

- ▶ Redesigned to mirror the contours of the new ergonomic seat plate, users of the soft seat will not lose any of the positioning benefits the seat plate offers.
- ▶ When the commode function is not required the soft seat insert can be fitted for additional comfort.
- ▶ Also available with a smaller aperture, the small soft seat caters for smaller or thinner adults or children.
- ▶ The Universal soft seat gives an alternative to the soft seat and soft seat insert when the commode function is not required or the user may have existing skin issues or skin breakdown, which would prevent the use of the insert.



▶ i) Standard soft seat




▶ ii) Standard soft seat Insert



▶ iii) Small soft seat



▶ iv) Universal soft seat

**i****j**

Seat and backrest accessories

Ergonomic soft seat

For users with existing pressure ulcers or skin breakdown.

- ▶ Manufactured from Viscoelastic (memory foam), allowing greater depth of immersion and a much softer feel.
- ▶ Optimal stability, comfort and positioning of the user due to the increased depth of immersion.
- ▶ As with all the soft seats, if the decision is made to use the ergonomic soft seat, the seat to ground height will increase. This may have an impact on those who are normally able to transfer themselves independently.
- ▶ Seat to ground height, footrest height and armrest height should all be checked when combining the existing seat plate with any of the soft seat accessories.

Seat and backrest accessories

Variable soft seat

Optimal access for hygiene purposes or for those who transfer sideways.

- ▶ The aperture of the variable soft seat can be positioned in four different directions.
- ▶ To fit the variable soft seat, the ergonomic seat plate must first be removed.
- ▶ Brackets for the variable seat are easily attached to the seat frame and secured with clips, after which the seat plate is clipped back into the desired position.

Seat and backrest accessories

Plastic backsupport and soft backrest cushion


The plastic backrest offers an alternative to the standard tension adjustable backrest when firmer support is required to achieve a more upright posture.

- ▶ The plastic backrest is not available or recommended on models with tilt function.
- ▶ The soft backrest cushion is designed for comfort, but can also help reduce the overall breadth or seat depth of the chair for smaller users or children.
- ▶ When using with thinner, taller users, adjust the tension adjustable backrest straps to ensure no loss of seat depth.



More Aquatec Ocean Ergo resources available...

Literature



1630242



1631734



1630142

Videos



YouTube



Visit www.invacare.co.uk to access all of the resources and much more!

References

- > **Jones, M. and Gray, S.** Assistive Technology: Positioning and Mobility. Meeting the Physical Therapy Needs of Children. FA Davis Company. 2005.
- > **Lustig, M., Levy, A., Kopplin, K., Ovadia-Blechman, Z. and Gefen, A.** Beware of the toilet: The risk for a deep tissue injury during toilet. Sitting. Journal of Tissue viability. 2017.
- > **Neville, L. Commissioned by James Leckey Design Limited.** The Fundamental Principles of Seating and Positioning in Children and Young People with Physical Disabilities. 2005.
- > **Purves, D., Augustine, G.J., Fitzpatrick, D., Katz, L.C., LaMantia, A., O McNamara, J. and Williams, M.** Neuroscience, 2nd edition. Sunderland (MA): Sinauer Associates. 2001.
- > **Romanelli, M., Clark, M., Gefen, A. and Ciprandi, G.** Science and Practice of Pressure Ulcer Management. Springer. 2018
- > **Wounds International.** International Review Pressure Ulcer Prevention, Pressure, Shear, Friction and Microclimate in Context. A Consensus Document. 2010.

Website References

- > <https://adrenalfatiguesolution.com/fight-or-flight-vs-rest-and-digest/>
- > http://www.dantest.com/dtr_ans_overview.htm
- > https://www.diffen.com/difference/Parasympathetic_nervous_system_vs_Sympathetic_nervous_system
- > <https://mobilitymgmt.com/articles/2017/09/01/seating-strategies.aspx>
- > <https://www.sciencedirect.com/book/9780323039307/primary-care-geriatrics>
- > <https://www.sis.se/api/document/preview/918591/ fight-or-flight-vs-rest-and-digest/>
- > www.spinalhub.com.au/what-is-a...bowel.../nervous-system-control-of-your-bowel
- > <https://www.squattypotty.com/blog/5-problems-with-sitting-on-your-toilet/>

Biography

Rhona Moot is a qualified Occupational therapist with 15 years' experience.

- She is currently Clinical Specialist for Invacare EMEA and has previously worked as Clinical Education Specialist for Seating and Positioning, Invacare EMEA. Prior to working at Invacare, Rhona worked for NHS Greater Glasgow and Clyde in Glasgow where she specialised in Wheelchair Assessment and Provision for Adults aged 16 years and older, whom has a range of physical disabilities.

Invacare Limited
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail: ordersuk@invacare.com
www.invacare.co.uk

All rights reserved. All information quoted is
believed to be correct at time of print.

Ocean Ergo Family
Clinical Introduction - UK - 02/2020
Available in Digital Format only.



Yes, you can.[®]

Safe
Patient
Handling