Invacare® Jasmine™

User's Manual	EN
Manuel d'utilisation	FR
Manual del usuario	ES
Gebrauchsanweisung	DE
Brukervejledning	NO
Manual do utilizador	PT
Gebruiksaanwijzing	NL

*This manual MUST be given to the user of the product.
BEFORE using this product, read this manual and save
for future reference.*

Yes, you can.®

User's Manual

EN

Stamp of the Distributor

Table of Contents

1.	General	4
	1.1. Symbols	5
2.	Safety	6
3	Technical data	7
4.	Components and functions	10
5.	Set up and installation	11
6.	Using the hoist	12
	6.1. Emergency, lowering, lifting and stop	15
7.	Maintenance and cleaning	16
8.	Troubleshooting	18
9.	Product disposal and recycling	19

Invacare®'s *Jasmine*™

Thank you for choosing *Invacare*®'s *Jasmine*™ mobile hoist. The *Invacare*® *Jasmine*™ is especially designed to lift the patient to and from wheelchairs, beds, toilets and floors. The *Jasmine*™ hoist offers easy handling and increased comfort within the nursing sector. The design of *Jasmine*™ provides a longer reach, wider leg span and a wider lifting range.

The electrical leg span combined with the large range of spreader bars and slings available from *Invacare*® makes it easy to adjust the hoist for optimum patient care.

Invacare® is certified according to DS/EN ISO 9001 and ISO 13485 which ensures that our customers are always supplied with products of uniform quality. Please see *Invacare*'s bulletin PU001 on the use of other manufacturer's slings with *Invacare* hoists.

Throughout the entire production process, our materials and components are quality controlled by the operators. A final test is made when the product is fully assembled. If the product does not correspond to the quality demands of *Invacare*®, it will not pass quality control.

In the event of a problem in connection with the delivered product, please contact your local *Invacare*® supplier.

Invacare® will take no responsibility if the product is used or assembled in any way other than stated in this user manual. Only accessories mentioned in this user manual may be used on the *Jasmine*™ mobile hoists. This product must only be operated by qualified staff who have received the necessary instruction and training.

Please read this user manual carefully before using the mobile hoist.

I. General

- *Jasmine*™ is CE-marked in accordance with directive 93/42/EEC concerning medical devices.
- *Jasmine*™ has been tested and approved according to EN/ISO 10535.
- Control unit and motors have been approved according to EN 60601.
- *Jasmine*™ has undergone a risk analysis according to EN/ISO 14971.
- Control unit, motors and hand control are IPX4 protected.
- Lifting capacity: 200 kg (31,5 stone).
- The hoist has an expected lifetime of 8 years.
- The hoist will not operate when in charge mode.
- Warranty - Frame 2 years, Actuator/ Motor 5 years, Battery 1 year.

Disconnect the charger cable before moving or using the hoist. Make sure the mains cable is not squeezed or damaged when assembling or dismantling the hoist or when moving the hoist.

If the functions of the hoist change, please see the section on "Maintenance and cleaning". This product must be serviced and maintained by qualified personnel.

If the hoist is used in the area of a shower or bath, ensure that the hoist is wiped clean of any moisture after use. **DO NOT** store the hoist in a damp area or in a damp condition. Periodically inspect all components of the hoist for signs of corrosion. Replace all parts that are corroded or damaged.

1.1. Symbols

The patient is not separated from the ground and the chassis:

Direct current:

Alternating current:

Double insulated:

200 kg (31,5 stone) max. load (SWL) = (Patient + sling + spreader bar):

The product should be recycled where possible:

Refer to user's manual:

2. Safety

Please be aware of the following possible risks when using the hoist:

The hoist must be carefully supervised if the patient is a child.

Please note that it is more difficult to manoeuvre the hoist on carpeted or uneven/textured surfaces as opposed to wooden or flat surfaces; this is due to the starting forces required to initiate the movement that may create additional friction exerted on the castors.

It is recommended not to move the user with the jib at maximum height. Take great care when moving the hoist on carpeted, wet, slippery, rough or uneven surfaces. Never use on slopes.

Caution! Never attempt to lift a patient outside the base area of the hoist. When lifting a patient from the floor, the patient must always be located between the two hoist legs, and never outside this area. To avoid injury, utmost care must be taken when raising and lowering a patient.

Invacare[®] accepts no liability for any use, change or assembly of the *Jasmine*[™] hoist other than as stated in this user manual.

There is a risk of squeezing and entrapment during operation and transport of the hoist. Never exceed the maximum recommended lifting capacity.

It is important that the cables do not become entrapped or jammed in moving parts and castors during operation and transport.


Important! *Invacare*[®] recommend that operation of any mobile lifting equipment is carried out by personnel with prior professional training in moving and handling techniques. *Invacare*[®] can provide Moving and Handling training in relation to *Invacare*'s hoist range. Please contact *Invacare* on 01656 776222 for further information.

Always choose the sling design and size according to the patient's weight, size and physical ability. Always ensure, that the sling is applied correctly. Be sure to check the sling attachments each time the sling is removed and replaced, to ensure that it is properly attached before the patient is removed from a stationary object (bed, chair or commode). Never leave a patient in the sling unattended. Please ensure that the hoist is not exposed to water, as this can cause damage.

3. Technical data


General specifications	
Highest position (max. height of CSP*) (k)	212 cm
Lowest position (min. height of CSP*) (l)	74 cm
Hoisting range (m)	138 cm
Lifting area (height range)	54 - 192 cm
Height to upper edge of legs (t)	12 cm
Min. free height (s)	4,5 cm
Minimum internal width (r)	53 cm
Internal width at maximum reach (q)	103 cm
Maximum Length of the base (u)	134 cm
Maximum Internal Length of the Base (v)	90.5 cm
Maximum External Width of Legs open (y)	128 cm
Total width (open) centre to centre of castors (p)	122.5 cm
Minimum Internal Width of Legs (open) (z)	118 cm
Maximum Internal Width of Legs (closed) (w)	69 cm
Minimum Internal Width of Legs (closed) (x)	56.5 cm
Castor size (FRONT/REAR)	10 cm
Sling Material	Polyester
Maximum lifting capacity (patient+sling+spreader bar)	200 kg (31,5 stone)
Total weight (weight out of carton)	43 kg
Weight, Mast including battery, spreader	21.2 kg
Weight, Leg section	21.8 kg
Maximum reach at 60 cm (a)	60 cm
Maximum reach from base (b)	66 cm
Reach from base with legs spread to 70 cm (c)	24.5 cm
Turning radius	153 cm

* CSP : Central Suspension Point

General specifications	
Operating forces of buttons Max.	5N
Battery (voltage output)	24V DC === max. 240VA
Charger Input (voltage supply)	100-240V AC ~ 50/60 Hz
Charger Output/Charging Time	29.5V DC Max 6 hours
Audio/Visual Low Battery Alarm	Yes
Motor Safety Devices	Anti-Entrapment
Working ability	*40 full lifts without battery charge with batteries at 50% of full capacity
Velocity of lifting and lowering	<0.15 m/s under maximum load and <0.25 m/s unloaded
Max. current input	Max. 400 mA
Operating temperature	5°C to 40° C
Air humidity	20% to 90% @ 30 not condensing
Sound pressure	< 50 dB
Atmospheric pressure	700 to 1060 hPa
Protection class, control unit	IPX4
Protection class, hand control	IPX4
Protection class, motor	IPX4
Insulation class	II - type B
Intermittens	10%, max, 2 minutes/18 minutes
Battery capacity	2,9 Ah
Manual emergency lowering	Yes
Electric emergency lowering/lifting	Yes/Yes
Expected lifetime	8 years

*NOTE: Varies depending upon load and stroke.

4. Components and functions

Accessories

It is compulsory to use original *Invacare*® Accessories which you can purchase through any *Invacare*® dealer.

For repair, please contact your local *Invacare*® dealer.

Service and repair must be carried out by a properly trained technician.

5. Set up and installation

Unpacking


1. Should the packaging be damaged upon receipt, each part of the hoist must be examined for visible defects. In the event of any damage, please contact *Invacare*[®].
2. Carefully take out and identify all the parts in the box.
3. The box contains the following parts:
 - 1 hoist
 - 1 user's manual
 - 1 battery
 - 1 mains cable
 - 1 spreader bar

To prevent self-discharging of the battery or accidental operation of the hoist, the emergency stop is activated during shipping. Please remember to deactivate the emergency stop and charge the battery prior to first use.

Assembling the hoist

The hoist in the *Invacare*[®] range are designed for easy assembly. The only tools required are a 1/2" spanner and a 3/16" Allen key.

1. Put the base on a level surface so that all castors make contact with the floor.
2. Lock the rear castors.
3. Loosen the bottom knob, but leave it screwed into the base.
4. Unscrew the top knob and remove it from the base.
5. Connect the mast assembly cable (not shown) to the base cable (not shown).
6. Pull the tube of the mast out of the hole in the base.
7. While supporting the mast assembly, tighten the bottom plastic knob.
8. Screw the top plastic knob into the hole in the base.
9. Using a 1/2" spanner and a 3/16" Allen key, mount the spreader bar to the moving part of the mast. Do not over-tighten the fixing nut.

Dismantling the hoist

The hoists are easy to assemble without the use of tools.

1. Put the base on a level surface so that all castors make contact with the floor.
2. Lock the rear castors.
3. Loosen the bottom knob, but leave it screwed into the base.
4. Unscrew the top knob and remove it from the base.
5. Pull the tube of the mast out of the hole in the base.
6. Disconnect the mast assembly cable (not shown) to the base cable (not shown).


The plastic bushing in the base that surrounds the mast must not be removed as this will cause instability of the mast. Activate the emergency stop before assembly or disassembly to prevent entrapment/squeezing. There is a possible risk of squeezing of both limbs and wires during assembly and disassembly of the hoist. Take the utmost care when lifting components during assembly - some parts are heavy. Always remember to adopt the correct lifting position.

Make sure there is an audible click when mounting the battery onto the (optional) battery charger to confirm proper mounting. Otherwise, injury or damage may occur.

Attaching the battery charger to a wall (applicable only to version with separate battery charger)

1. Place the battery charger mounting bracket on the wall at the desired position.
2. With a pencil, mark the middle hole position.
3. Measure down 165 mm from the pencil mark and drill one mounting hole.
4. Install the bottom mounting screw until there is an approximate 3 mm gap between the screw head and the wall.
5. Put the battery charger mounting bracket onto the bottom mounting screw.
6. Drill the other two mounting holes.
7. Screw the mounting screws through the battery charger mounting bracket and into the wall. Tighten securely.
8. Put the battery charger in place on the mounting bracket.
9. Plug the battery charger into an electrical outlet.


6. Using the hoist

Lifting capacity: 200 kg (31,5 stone).
(patient+sling+spreader bar)

Using the hand control

1. Press “arrow up” on the hand control to raise the jib.
2. Press “arrow down” on the hand control to lower the jib.
3. Operation of the electric legs.

NB - Operating forces of buttons: Max. 5 N.


WARNING! When using the Jasmine™ in conjunction with beds or wheelchairs, be aware of the position of the lift in relationship to those other devices so that the lift does not become entangled.

LED diodes on the control unit

The left green diode will light continuously when the control unit is connected to the mains. The yellow right diode will light continuously during charging, and go out when the battery is fully charged. There is no indication when hand control buttons are pressed (except from the remaining battery capacity display or when the emergency lowering/raising system is activated).

Brakes

1. Stand behind the hoist and take a good hold of the handle bar.
2. Press the pedals downwards with the foot on the back castors to activate the brake.
3. Press the pedals upwards with the foot to release the brake.

Never try to activate/deactivate the brakes by hand - risk of squeezing !

Moving the hoist

1. Stand behind the hoist and take a good hold of the handle bar.
2. Release the brakes of the back castors.
3. Now the hoist can be pushed or pulled to the desired location.

Note: To negotiate raised obstacles during movement of the unloaded hoist, it is recommended to pull the hoist backwards. This operation should only be performed when the hoist is in the unloaded condition.

Use handle bar on the mast at all times to push or pull the patient hoist.

Turning in place

The hoist can be turned (rotated) in place for transfers with limited floor space.

Charging the battery

The hoists are equipped with an internal charger. It is recommended to charge the batteries regularly to ensure optimal use of the hoist and prolong the the batteries. Furthermore, it is recommended to charge the batteries before first use.

The control unit is equipped with a sound signal, which will beep when operating with low battery capacity. It is recommended to charge the batteries as soon as the sound signal is heard. In addition, *Jasmine™* is equipped with a display, indicating the approximate remaining battery capacity, when pressing any hand control button. When installing the battery onto the control unit, it takes a few seconds before the hoist is operational.

How to charge the batteries:

1. Connect the mains cable to the control unit and plug it in. It takes about 4 hours to charge the batteries. Charging must take place in a room with good air ventilation. The charger stops automatically when the batteries are fully charged.
2. Remember to disconnect the charger cable before using the hoist again.

The emergency stop must not be activated - otherwise it will be impossible to charge the batteries. While charging takes place the hoist cannot be used. Do not use or move the hoist without unplugging from the socket outlet. Do not attempt to use the hoist if the battery housing is damaged. Replace a damaged battery housing before further use.


Charging the Battery with battery charger

1. Lift UP on the handle on the back of the battery.
2. Lift the battery up and away from the control box.
3. Place the battery on the battery charger. Push the top of the battery against the mounting bracket until there is an audible click.

NOTE: The charge LED will illuminate. When charged, the LED will stop illuminating.

NOTE: It will take approximately four hours to charge a battery that requires a full charge.

4. Lift UP on the handle on the back of the battery.
5. Lift the battery up and away from the battery charger.
6. Place the battery on the control unit. Push the top of the battery against the mounting bracket until there is an audible click

CAUTION: Make sure there is an audible click when mounting battery on the battery charger. Otherwise, injury or damage may occur.

Fitting the sling

Slings designed for 2- or 4-point spreader bars can be used with this mobile hoist, please refer to *Invacare*®'s sling brochure for details. For sling compatibility, please see *Invacare*®'s bulletin PU001. Copies can be obtained by calling Customer Services on 01656 776222.

Note: The size of the sling and the width of the spreader bar should correspond.

The *Invacare*® Universal Sling can be used for highly physically dependent patients; however, they must have a good degree of head control. The sling allows for easy access for personal hygiene and it is easily applied in both sitting and lying position.

1. Select the sling best suited for the patient in question. Do not attempt to lift any user without prior professional training. If in doubt ask your *Invacare*® dealer for advice.
2. Fit the sling to the patient as instructed in the *Invacare*® sling user's guide. Take great care to ensure that the sling is fitted correctly, so that the patient is comfortable and safe.
3. Manoeuvre the hoist, so that the spreader bar can be lowered to the sling loops as close as possible.
4. Widen the chassis legs for increased stability.
5. Lower the spreader bar; attach the 4 loops of the sling by sliding the loops into the hooks.
6. If the patient is in a sitting position, lower the spreader bar to patient's chest, then move the hoist nearer to the patient. Ensure that brakes are not engaged prior to lifting or lowering.
7. Sitting in chair or wheelchair Lean the client forward ensuring that they are well supported to reduce risk of falling and to promote feelings of security. Slide sling down behind the back of the patient in a central position with the aperture apex aligned with the base of the spine. Sit the patient back on against the sling. Pass the leg support towards the front of the chair to each side of the thigh and under the thigh.
8. Lying in bed, on floor or in bath Lean the patient forward ensuring that they are well supported to reduce risk of falling and to promote feelings of security. Slide sling down behind the back of the patient in a central position with the aperture apex aligned with the base of the spine. Sit the patient back on against the sling.


Caution! Do not lift or lower a patient with the brakes applied. Always let the hoist find the correct centre of gravity. Never attempt to lift a patient outside the base area of the hoist. When lifting a patient from the floor, the patient must always be located between the two hoist legs, and never outside this area. To avoid injury, utmost care must be taken when raising and lowering a patient. **Never** exceed the recommended maximum hoist weight and never use the hoist in wet areas.

6.1. Emergency, lowering, lifting and stop

Electrical emergency lowering

If the hand control fails, the jib can be lowered by using the circular recessed button for emergency lowering. This is located at the front of the control unit.

The jib will be lowered as long as the button is pressed, and should be operated by means of inserting a pointed object, such as a pencil into the circular recess.

Electrical emergency lifting

If the hand control fails, the jib can be raised by using the button for emergency lifting. This is located at the front of the control unit, and should be operated with a pointed object, such as a pencil. The jib will be lifted as long as the button is pressed.

Manual emergency lowering

In case of partial or total power failure, or if the battery runs down while using the hoist, *Jasmine*TM is equipped with a manual emergency lowering system located at the bottom of the motor. The weight is preset for 75 kg. If the patient's weight is more or less than 75 kg, it might be necessary to adjust the manual emergency lowering.

Adjustment of the manual emergency lowering

If the lowering speed has to be adjusted according to the patient's weight, then the screw in the red release handle must be adjusted. Loosen to increase the speed and tighten to decrease the speed.

Note: Manual emergency lowering is only possible when a patient is sitting in the hoist.

Emergency stop

In case of emergency: Press the red button for emergency stop. The emergency stop can be reset by turning the stop button

Note! The mechanical emergency lowering system will only operate when there is a patient in the hoist.
The mechanical emergency lowering system has to be adjusted according to the patient's weight.

7. Maintenance and cleaning

At normal daily operation, a service check-up should take place every year, according to the maintenance chart.

When performing annual or regular maintenance, all parts designed to carry load must be as a minimum tested with maximum load. All safety features must be checked according to EN ISO 10535: 2006 Annex B.

LOLER statement.

The UK Health and Safety Executive's Lifting Operations and Lifting Equipment Regulations 1998, require any equipment that is used in the workplace to lift a load be subject to safety inspection on a six monthly basis. Please refer to the HSE web site for guidance www.hse.gov.uk

It is the responsibility of the person responsible for the equipment to ensure that LOLER regulations are adhered to.

The hoist must be stored at normal room temperature. If it is stored in a damp, cold or wet environment then the motor and other mounting parts may be prone to corrosion.

- The hoist has been designed so that it requires a minimum amount of maintenance.
- It is recommended that the batteries are charged frequently, if possible every night, as it will extend the lifetime of the batteries and will ensure high performance.
- It is recommended that the hoist is cleaned after use.
- Please check that the spreader bar connections, control unit and motors are in good working order and not damaged in any way. If there is any damage to the hoist, refrain from using it, and contact your local *Invacare*® dealer before further use.

Wear and damage of loading parts

It is necessary to maintain and check all parts that are exposed to static or dynamic strain, such as the sling, the spreader bar and the turning points, do not have fractures, are fragile, askew or damaged.

Damaged parts must be replaced before further use.

Cleaning

- To prevent cross-infection, the hoist must be cleaned after use.
- The hoist must be wiped with a moist firmly wrung cloth with ordinary household disinfectants.
- Only use officially approved disinfection detergents.
- Dry the hoist carefully after cleaning.
- Never use acids, alkaline or solvents for cleaning the hoist.

Motors, control unit and mounting parts can be destroyed if the hoist is cleaned in any other way than stated above.

Maintenance chart

Only personnel having received the necessary instruction or training by Invacare® must perform service and maintenance on Jasmine™.

Serial number (located on the hoist base): _____

Date:	Initials:				
Visual inspection of all hoist parts					
Inspection of the spreader bar - weldings and assemblies					
Inspection of the motors - undamaged					
Inspection of cables and plugs - undamaged					
Inspection of emergency stop, electrical emergency lowering and raising					
Inspection of the castors - attachment, braking ability and free rolling					

Invacare® offers courses in service and maintenance of the hoist.


For all servicing and maintenance requirements, please contact your local Invacare® dealer.

IMPORTANT: Please see LOLER statement. Page 16, section 7.

JUMBO Care functions

For a safe maintenance, the option JUMBO Care offers five (5) main functions:

- 1) Service intervals
- 2) Service information on PC
- 3) Actuator life time
- 4) Resetting of service interval
- 5) Actuator exchange

Checking the Service Light

Note: Each time the hoist is assembled, and before using the hoist, the service light should be checked.

Initial hoist assembly or assembly within 12 months of use

1. Turn on the hoist.
2. Examine the control box to see if the service light is flashing.
3. Perform one of the following:
 - A. Service Light is Flashing - The service light needs to be reset. Locate the handset. Press and hold the UP and DOWN buttons at the same time for five seconds (1-2).

Note: You will hear a sound when the service light has been reset.

B. Service Light is Not Flashing - The hoist is ready for use.

Charge the battery if necessary. Assembly after more than 12 months of use

1. Turn on the hoist.

2. Examine the control box.

3. Perform one of the following:

A. Service Light is flashing - The hoist requires service.

Contact your local Invacare dealer for service.

B. Service Light is not flashing - The hoist is ready for use.

Charge the battery if necessary.

8. Troubleshooting

Only personnel having received the necessary instruction or training by *Invacare*® must perform service and maintenance on *Jasmine*™.

Symptom	Possible cause	Remedy
The mast seems loose	The mast screws has become loose	Tighten the screw
Castors are noisy	Dust or dirt in the castors	Clean the castors
Hinge between jib and mast is noisy during operation	Missing lubrication	Lubricate the hinge
The motor is not running	Hand control or motor is not connected	Insert the plugs fully
	No power on the battery	Charge or replace the battery
	The red emergency stop is activated	Turn the red button clockwise to deactivate
Motor noise, but no movement in the piston rod	Motor is damaged	Replace the motor
The piston rod only moves in and not out	Motor is damaged	Replace the motor
The jib motor stops during a patient transfer	Max. load is exceeded	Reduce the load (and the hoist will function normally)

Contact your dealer if the above does not solve your problems.

Lubricate with medically clean oil, e.g. Kemitura Kem Lub KEM-WO 50, order no. 813239.

9. Product disposal and recycling

This product has been supplied from an environmentally aware manufacturer that complies with the Waste Electrical and Electronic Equipment (WEEE) Directive 2002/96/CE.

This product may contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to legislation.

The 'crossed out wheelie bin' symbol is placed on this product to encourage you to recycle wherever possible.

Please be environmentally responsible and recycle this product through your recycling facility at the end of its life.

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INVACARE AG
Invacare[®] Rehabilitation Equipment (Suzhou) Co., Ltd