Invacare® Matrix® accessory

en	Headrest pads / hardware User Manual	1
de	Kopfstützenpolster/Befestigungsteile Gebrauchsanweisung	4
fr	Housses de rembourrage/matériel des appui-têtes Manuel d'utilisation	6
it	Cuscini del poggiatesta / hardware Manuale d'uso	8
es	Piezas/almoHADILLAS del reposacabezas Manual del usuario	9
pt	Almofadas de apoio de cabeça/ferragens Manual de utilização	12
nl	Hoofdsteunkussens/hardware Gebruiksaanwijzing	14
da	Puder/fastgørelsesanordning til hovedstøtte Brugsanvisning	16
no	Hodestøtteputer/maskinvare Bruksanvisning	18
sv	HuvudstödSDynor/monteringsenheter Bruksanvisning	19

This manual MUST be given to the user of the product.
BEFORE using this product, read this manual and save for future reference.

Yes, you can.®

I General

I.1 Images


Fig. 1


Fig. 2

1.2 Symbols in this user manual

In this User Manual warnings are indicated by symbols. The warning symbols are accompanied by a heading that indicates the severity of the danger.

WARNING

Indicates a potentially hazardous situation which if not avoided could result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation which if not avoided could result in product damage, minor injury or both.

IMPORTANT

Indicates a hazardous situation which if not avoided could result in damage to the product.

Gives useful tips, recommendations and information for efficient, trouble-free use.

This product complies with the directive 93/42/EEC for medical products. The launch date for this product is specified in the CE declaration of conformity.

Manufacturer

1.3 Warranty

We provide a manufacturer's warranty for the product in accordance with our General Terms and Conditions of Business in the respective countries. Warranty claims can only be made through the provider from whom the appliance was obtained.

1.4 Standards and regulations

Quality is of the utmost importance for our company; all processes are based on the ISO 9001 and ISO 13485 standards.

This product features the CE mark, in compliance with the Medical Device Directive 93/42/EEC Class I.

The product meets the requirements of the standards DIN EN 12182 (Technical aids for disabled persons - general requirements and test methods), DIN EN ISO 14971 (Medical devices - Application of risk management to medical devices) and DIN EN ISO 10993-1 (Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management system).

Invacare® is continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum. We only use REACH compliant materials and components.

For further information please contact Invacare® in your country (for addresses, see back page of this manual).

1.5 Intended use

The Matrix® Elan headrest hardware and headrest pads are designed to offer comfort and posterior head support. It will provide simple, comfortable and effective support helping to maximise function and seating tolerance.

1.6 Scope of delivery

The following components are included in the scope of delivery:

- Matrix® Elan Hardware (MEHW, post length 356 mm) or Matrix® Elan Mini Hardware (MEMH, post length 254 mm) (see Fig. 1)

- Headrest pad (see Fig. 2 and following table for options)

Option	Designation
(A)	Elan Standard Pad 150 mm x 90 mm – SMALL
(B)	Elan Standard Pad 250 mm x 130 mm – MEDIUM
(C)	Elan Standard Pad 360 mm x 130 mm – LARGE
(D)	Elan Occipital Pad 230 mm x 180 mm
(E)	Elan Occipital Pad 300 mm x 200 mm
(F)	Elan 4-Point Pad 280 mm x 260 mm

- User manual

1.7 Service life

The expected service life of this product is three years when used in accordance with safety instructions and intended use, stated in this manual. The effective service life can vary depending on frequency and intensity of use.

2 Setup and usage

2.1 Installing the mounting hardware / headrest pads (Fig. 1)

WARNING!

Risk of serious injury

Injury may occur during use of the wheelchair if a headrest is wrongly adjusted or not installed correctly. The users neck may become hyperextended in the event of a collision, or when the system is tilted or reclined.

- This item must be fitted by an experienced clinician or technician.
- Before each use, check that all screws are securely fastened.

WARNING!

Risk of severe injury when transporting the client within a wheelchair fitted with an Invacare® Matrix® ELAN headrest in a motor vehicle.

- Wherever practicable Invacare recommends that the client is transferred out of the wheelchair into a vehicle seat using a lap and diagonal seat belt or an appropriate child safety seat.
- Invacare® Matrix® ELAN headrests may only be used in conjunction with wheelchairs that comply with the performance requirements of RESNA WC19 or ISO 7176-19.
- For safe transportation of a client within a wheelchair in a motor vehicle, refer to the user manuals of the wheelchair and the tie-down and occupant restraint system (WTORS).

1. Remove the headrest assembly and components from the packaging.

WARNING!

Risk of injury

- Inspect the headrest hardware to ensure that there are no missing or damaged parts.

WARNING!

Risk of injury

- Before drilling or modifying any back shell, always contact the back manufacturer to ensure the back is recommended for a headrest application.

Tools:	Allen keys 2.5 mm / 4 mm / 5 mm (provided) Hacksaw, Power drill, 6 mm drill bit
--------	--

- Secure the headrest clamp **B** onto the back shell using the hardware provided. Install the clamp using the existing mounting holes in the back shell (or adapter plate). If mounting holes are not provided, holes may be carefully drilled in the desired location using the mounting plate as a template.

 An optional spacer plate **D** with extended mounting screws may be installed for additional clearance when needed.

- Install the desired headrest pad **I** onto the MatrX® Elan hardware using the headrest mounting ring **F**, upper multi-angle rotational pivot **G** and the 3 screws **J** provided.

 Tighten the hardware just enough to secure the headrest pad in position, but still allow for final (fine tune) angle adjustments via the upper pivot mechanism.

- Loosen/remove the lower D-ring **C** from the vertical mounting post **A** using the 2.5 mm Allen key provided.
- Slide the vertical mounting post into the clamp assembly **B**.
- Adjust the overall height of the headrest pad to the desired position by positioning the upper D-ring **K** on the vertical mounting post as required.
- Tighten knob **L** to secure the MatrX® Elan hardware in place.
- Adapt the headrest to the desired depth and angle according to the needs of the user by adjusting the 3 linkages **H** of the dual link adjustment assembly using the 4 mm Allen key, in combination with the upper and lower multi-angle (rotational) pivots **G** and **E** using the 5 mm Allen key.

! IMPORTANT
Do not fully remove the dual link adjustment hardware.

– For angle adjustment loosen the screws in the linkages **H** for only approx. 1/2 turn. Optional lever handles **N** are available for easier adjustment of the upper and lower linkages. For installation carefully remove only the main screw from the upper or lower linkage (leaving the hex insert in place), then screw the lever into the hex insert.


- Once the final position is established, fully tighten all linkages and headrest hardware.
- Secure the Lower D-Ring onto the vertical post so that it rests flush with the bottom of the clamp assembly.

 If necessary, the overall length of the vertical post can be modified (cut) in order to prevent any interference along the back of the seating system. After cutting the post, the plastic end cap **M** in the bottom of the post can be removed and transferred to the newly cut post.

3 Cleaning and care

! IMPORTANT!
All cleaning agents and disinfectants used must be effective, compatible with one another and must protect the materials they are used to clean.
– For further information on decontamination in health care environments, please refer to 'The National Institute for Clinical Excellence' guidelines on Infection Control www.nice.org.uk/CGI39 and your local infection control policy.

- Clean the pads with a neutral detergent and hot water and / or a single use wipe.
- Clean the mounting hardware with common household cleaners and a soft cloth.
- Store in a dry environment.

 Should the product become contaminated, ask your infection control officer for advice. Alternatively, contact Invacare for their information on the decontamination protocol.

4 Materials

Assembly hardware	Aluminum / Steel
Headrest pad cover	100% Interlock Polyester (black) with Polyurethane Coating
Headrest pad foam	Moulded High Resilience (HR) Polyurethane Foam (green)

All materials used are corrosion resistant and not made with natural rubber latex.

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Advena Ltd.