



Invacare® EC-Track™, Gantry, Robin® EVO

en Periodic Inspection Checklist



PROVIDER: Keep this manual.
The procedures in this manual **MUST** be performed by a
qualified technician.



Yes, you can.®

Periodic Inspection Checklist | EC-Track System, Gantry and Robin EVO

In accordance with ISO 10535.



CAUTION!

Risk of injury and damage to property

Failure to ensure an adequate **Periodic Inspection** of the product can compromise the safety of people and property.

- A periodic inspection of the product must be performed at least every **12 months** unless otherwise stated in local requirements.
- The frequency of inspection must be increased if the product is continuously exposed to high humidity, high condensation and corrosives (e.g., chlorine or ammonia gases), to prevent product impairment.
- Inspections **MUST** be performed by a qualified person who is acquainted with the design, use and care of the product.
- The inspection record **MUST** be recorded in the logbook and reported to the owner.

Is the EC-Track on the installation?

Yes (Complete the following section):

The Invacare EC-Track System installation label shows the following information:			
Maximum SWL		Installation Date	
Label Number (LN)			

No (Proceed to **2. Checkpoints | Robin EVO**).


1. Checkpoints EC-Track and Gantry	✓	✗	Comments
1.1 The installation label is present, visible and readable.			
1.2 All connection points, brackets to structure are tightened with recommended torque. Retighten, if necessary (if NOT accessible please see point 1.21 on load test).			
1.3 All connection points, brackets to rail system are tightened with recommended torque. If the system is the Gantry, all the connections are tightened with recommended torque. Retighten, if necessary (if NOT accessible please see point 1.21 on load test).			
1.4 All rail ends are secured with end stops, sufficiently tightened. Retighten, if necessary.			
1.5 End stops don't show any sign of damage, wear and/or deformation.			
1.6 Brackets don't show any sign of damage, wear and/or deformation. If the system is Gantry, connections don't show any sign of damage, wear and/or deformation.			
1.7 Trolleys don't show any sign of damage, wear and/or deformation.			
1.8 The trolley (and if applicable, traverse trolleys) run smoothly and silently in the entire track system.			
1.9 Check and adjust if required the trolley friction brake.			
1.10 End Covers are assembled (If applicable).			
1.11 Rail systems with track connections: The connections are aligned. Re-align, if needed.			

Transit Coupling / Track Switch (if applicable)					
1. Checkpoints EC-Track and Gantry	✓	✗	Comments		
1.12 Lock fixtures all OK (secure to tracking, each claw moves independently, correctly aligned with tracking slots, no obstructions to prevent correct function, safety plate is fitted).					
1.13 Check traverse trolley mounting fixtures.					
1.14 The tracking is aligned. Re-align, if needed.					
1.15 The gap between the track and the transit coupling/track switch is max. 3 mm.					
1.16 Function of transit coupling / lock fixture OK (10 times — disable left / right claw x 5).					
1.17 Function test performed under load through the track switch of the system (5 cycles with the maximum SWL).					
1.18 Inlay plate is fitted (L and M rail).					
Load test					
1.19 The track system is leveled and all fixings are sufficiently tightened before carrying out a load test.					
1.20 If all wall and ceiling brackets are accessible, and sufficiently secured a load test must be performed with the maximum SWL (200 kg). The load must be passed through the entire rail system.					
1.21 If all wall and ceiling brackets are NOT accessible a load test must be performed with 300 kg – 1,5 × maximum SWL (200 kg). This test must take a minimum of 20 minutes. The load must be suspended under each bracket/fixing point and passed through the entire rail system.					
1.22 The rail system remains the same height after the load test.					
1.23 The rail system remains leveled after the load test.					
1.24 Place a label with the date for next inspection (MM-YYYY) on the Tracks. Label Example:					
<table border="1" style="width: 100%;"> <tr> <td style="background-color: #cccccc;">Date for Next Inspection</td> </tr> <tr> <td style="height: 20px;"></td> </tr> </table>				Date for Next Inspection	
Date for Next Inspection					

Does this inspection include the Robin Evo?

- Yes (Complete the following section).
- No (Skip to the **3. Final Steps** and make sure to check the hoist’s manufacturer **Instructions / Installation Manual**).

2. Checkpoints Robin EVO	✓	✗	Comments
2.1 Check for broken, missing or misaligned parts.			
2.2 Product Label is present and legible (write SN on the comments section).			
2.3 The hoist is marked with the maximum SWL.			
2.4 The emergency stop and reset work correctly.			
2.5 The emergency lowering function and reset work correctly.			
2.6 The emergency lowering end stop switch operates correctly.			
2.7 The strap reset procedure works correctly.			
2.8 The full length of straps (tapes) have no signs of wear, fraying and/or loose stitching.			
2.9 Check if the end stop strap switch is working correctly (via hand control).			

2. Checkpoints Robin EVO	✓	✗	Comments		
2.10 Hooks and hook fixation (screws and pins) do NOT have signs of wear/damage that could endanger the lifting.					
2.11 Check suspension and cross pin.					
2.12 Steel and plastic gearwheels do NOT have signs of wear and damage. If gearwheels are dry, apply gearwheel grease. Refer to — step 5.					
2.13 Check tape guide cylinders.					
2.14 The hand control is connected correctly.					
2.15 The hand control is functional and in good conditions.					
Electrical Connections					
2.16 All cables are intact.					
2.17 No signs of rupture on cables and on isolation.					
2.18 Clamping and mechanical function of micro switches work correctly.					
2.19 The charger does not have signs of damage and is functional.					
2.20 Check that battery is charged when connecting the hand control to the charging station.					
Load Test					
2.21 Perform a full lifting cycle with the maximum SWL (200kg).  Note: Must be effortless and without jarring sound.					
2.23 Stick a label with the date for next inspection (MM-YYYY) on the Robin EVO hoist. Label Example: <table border="1" style="width: 100%;"><tr><td style="background-color: #cccccc;">Date for Next Inspection</td></tr><tr><td> </td></tr></table>	Date for Next Inspection				
Date for Next Inspection					

3. Final Steps	✓	✗	Comments
3.1 The owner has all the necessary User Manuals .			

Overall assessment

<input type="checkbox"/> Pass No safety or functional deficiencies were detected.	<input type="checkbox"/> Fail Safety or functional deficiencies were detected, and corrective actions are required. Comments:
--	---

Date of next inspection: _____

Service Company Name: _____

Inspection done by: _____

Date / Signature: _____

**United Kingdom & Ireland:**

Invacare Limited
Pencoed Technology Park, Pencoed
Bridgend CF35 5AQ
Tel: (44) (0) 1656 776 200
uk@invacare.com
www.invacare.co.uk

Belgium & Luxemburg:

Invacare nv
Autobaan 22
B-8210 Loppem
Tel: (32) (0)50 83 10 10
Fax: (32) (0)50 83 10 11
marketingbelgium@invacare.com
www.invacare.be

Danmark:

Invacare A/S
Sdr. Ringvej 37
DK-2605 Brøndby
Tel: (45) (0)36 90 00 00
Fax: (45) (0)36 90 00 01
denmark@invacare.com
www.invacare.dk

France:

Invacare Poirier SAS
Route de St Roch
F-37230 Fondettes
Tel: (33) (0)2 47 62 64 66
contactfr@invacare.com
www.invacare.fr

Italia:

Invacare Mecc San s.r.l.,
Via Marco Corner, 19
I-36016 Thiene (VI)
Tel: (39) 0445 38 00 59
servizioclienti@invacare.com
www.invacare.it

Nederland:

Invacare BV
Galvanistraat 14-3
NL-6716 AE Ede
Tel: (31) (0)318 695 757
nederland@invacare.com
www.invacare.nl

Norge:

Besøksadresse:
(Office addresses)
Invacare AS
Brynsveien 16
0667 Oslo
Tel: (47) 22 57 95 00
norway@invacare.com
www.invacare.no

Norge:

Vareleveringsadresse:
(Storage / Teknical dep)
Østensjøveien 19
0661 Oslo
eknisk@invacare.com
www.invacare.no

Sverige:

Invacare AB
Fagerstagatan 9
S-163 53 Spånga
Tel: (46) (0)8 761 70 90
Fax: (46) (0)8 761 81 08
sweden@invacare.com
www.invacare.se



Invacare UK Operations Limited
Unit 4, Pencoed Technology Park,
Pencoed
Bridgend CF35 5AQ
UK



Invacare Portugal, Unipessoal, Lda
Rua Estrada Velha 949
4465-784 Leça do Balio
Portugal



60147935-A 2024-09-24



Making Life's Experiences Possible®



Yes, you can.