



Simply mobile.

# scalamobil® scalacombi®

Mobile stair climbers



Using stairs  
the easy way.



*"Finally we can go anywhere with the wheelchair!"*

# Finally no more obstacles in the way – thanks to scalamobil and scalacombi

## Climbing stairs easily and safely – proven thousands of times

Every single step becomes an obstacle when an illness or handicap stops you climbing the stairs on your own: at home, at the surgery, in public buildings or when visiting friends and family. The mobile stair climbers from Alber – scalamobil and scalacombi – remove the obstacles from your path.

## And this is how it works...

Two electrically driven pairs of wheels continually climb one step after the next – upwards and downwards. This is ensured by a gentle yet reliable climbing mechanism – patented by Alber. When descending the stairs, four safety brakes stop scalamobil and scalacombi automatically at the edge of each step. The carer selects the direction and speed. When scaling the stairs, he or she only has to keep the stair climber in balance.

## The product benefits at a glance

- Also able to cope with narrow and spiral staircases – at home and elsewhere
- Functions on any floor covering – whether wood, stone, metal or carpeting
- No damage to stairs or floor covering
- ErgoBalance grips for safe operation with minimum exertion
- Climbs up to 300 steps with a single battery charge
- With integrated safety sensors



Go to the film here

You can also watch a film about scalamobil and scalacombi on the Internet at [www.alber.de](http://www.alber.de).








# scalamobil®


## Climb the steps in your own wheelchair

The scalamobil fits most common wheelchair models. The great advantage: you can stay seated in your wheelchair to climb the stairs.



①

A special **bracket** is fitted to your wheelchair and the scalamobil is attached when needed. This device does not restrict the folding mechanism or other wheelchair settings.



②

It takes just a few twists of the wrist to attach the **scalamobil** to the wheelchair bracket and secure it.



**Universal adapter scalaport** for the transport of various types of wheelchairs without any special brackets. Particularly suitable for transport services or public buildings.



*“Safely and comfortably up and down all types of stairs – this is quality of life”*

# scalacombi®

## Climbing the stairs with an integrated seat unit

Perhaps you need assistance to climb the stairs, but don't really need a wheelchair. The stair climber scalacombi will get you upstairs and downstairs with no effort at all. The compact dimensions and steerable front wheels make the scalacombi extremely versatile. It can be easily manoeuvred even in very small rooms and on narrow stairways.

When not in use, the scalacombi can be folded up to save space:

scalamobil and scalacombi also come as special versions for particularly high steps (up to 25 cm/9.8 inches). There is also a comprehensive range of accessories for all versions (e.g. belt systems, headrests).





# scalamobil and scalacombi – their particular strengths

## Intelligent and safe

The highest degree of safety and comfort when climbing stairs thanks to intelligent electronics: scalamobil and scalacombi can be programmed to adapt to the character of your stairs at home and will then climb them particularly gently and comfortably.

## Versatile

scalamobil and scalacombi adapt to your needs – not the other way around. The grip unit can be adjusted in height and width and thus adapts optimally to the person operating it and to the stairs.

## Easy to operate

With excellent ergonomics, all operating elements are in exactly the right place and the stair climber stays in an optimal position in the hands. This is where Alber's more than 25 years of experience in the development of mobile stair climbers can really be felt and experienced.



**ErgoBalance grips:** non-slip with a convenient shape



**Display showing** all functions at a glance



**Comfortable cushioning** for resting the operating unit on the thigh while climbing the stairs





*"We can now cope with any staircase even when we're out and about!"*

## Convenient technology on the go



scalamobil and scalacombi can be used everywhere. Your stair climber can be dismantled in a flash into three compact parts and will even fit in small car boots.

The battery pack can be removed with one hand, is less than 3 kg (6 lb) light and has enough energy for up to 300 steps with just one battery charge.

■ ■ ■ MADE IN GERMANY

Made in Germany by the market leader for add-on drive systems.



## scalamobil and scalacombi at a glance. Common questions and answers.

### **“Will the scalamobil also fit on my wheelchair?”**

Yes, the scalamobil fits on nearly all models of wheelchair by means of a small bracket that is mounted just once on the existing wheelchair.

### **“Do scalamobil and scalacombi also work on spiral staircases?”**

Yes, and even on different materials such as carpet, wood, stone or metal.

### **“Could my stairs be damaged?”**

No, as a result of the unique climbing mechanism and thanks to the abrasion-resistant tyres, the steps are optimally protected.

### **“How are scalamobil and scalacombi driven?”**

The quiet electric motors run on a powerful rechargeable battery that is kept operational using a charger.

### **“How far can I travel with a single charge of the battery?”**

As a rule, one charge of the battery is enough to climb 300 steps – for a person weighing approx. 70 kg (154 lb). That is equivalent to about 20 floors.

### **“Is a lot of strength needed to operate it?”**

At the push of a button, the electric motor takes care of all the climbing work. Guiding the scalamobil or scalacombi requires very little strength.

### **“What is the maximum permitted height and depth of a step?”**

The steps where a scalamobil is used should be at least 11 cm (4.3. inches) deep and not higher than 20 cm (7.9 inches). As an option we also offer a step height extension up to 25 cm (9.8 inches).

### **“How heavy can a person be who wishes to use the scalamobil or scalacombi as an aid to go up and down stairs?”**

The person's weight should not exceed 140 kg (308 lb)\*.

\* The maximum user weight may vary depending on the used wheelchair and the suitability of the person operating the device.

### **“Can the scalamobil and scalacombi also be used outside the house?”**

Basically, the stair climbers can be used on all kinds of steps, whether outdoors or in public buildings. The scalamobil can be dismantled into three parts; it is easy to transport and fits in every car boot.

### **“How much do scalamobil and scalacombi weigh?”**

scalamobil drive unit:	17.9 kg (39.43 lb)
scalacombi drive unit (incl. seat unit):	29,9 kg (65.9 lb)
Battery pack:	2.7 kg (5.9 lb)
Grip unit:	4.9 kg (10.8 lb)

scalamobil und scalacombi are approved medical device. Your medical supplies dealer will happily answer any questions regarding reimbursement.

scalamobil und scalacombi conform to Directive 93/42/EEC concerning medical devices and were voluntarily audited by TÜV Süd according to the current version of the product standard ISO 7176-28.

Our claim to quality is your guarantee of safety! Our company is certified to the latest quality management standard (including the respective latest amendments) set out in ISO 13485 for medical devices. With our environmental management system, certified to ISO 14001, we ensure the environmentally compatible manufacture of our products.

### **Weights scalamobil /scalacombi**

Max. user weight:	<b>140 kg (308 lb)*</b>
Max. total weight:	<b>180 kg (396 lb)*</b> <small>(scalamobil/scalacombi, wheelchair, occupant)</small>



**Alber GmbH**

Vor dem Weißen Stein 21  
72461 Albstadt, Germany  
Phone +49 07432 2006-0  
Fax +49 07432 2006-298  
info@alber.de

[www.alber.de](http://www.alber.de)