Izzo H340

en	Drop arm commode	
	User Manual	1
de	Nachtstuhl mit absenkbaeren Armen	
	Gebrauchsanweisung	4
fr	Chaise toilettes à bras rabattables	
	Manuel d'utilisation	7
it	Comoda con braccioli ribaltabili	
	Manuale d'uso	9
es	Inodoro con reposabrazos abatibles	
	Manual del usuario	11
pt	Sanita com braços amovíveis	
	Manual de utilização	13
nl	Toiletstoel met beweegbare armleuningen	
	Gebruiksaanwijzing	15

en

I General

I.1 Images


Fig. 1


Fig. 2

This manual **MUST** be given to the user of the product.
BEFORE using this product, read this manual and save for future reference.

Yes, you can.®


Fig. 3


Fig. 4


Fig. 5

Fig. 6


Fig. 7


Fig. 8

1.2 Symbols in this user manual

In this user manual, warnings are marked with symbols. The title next to the warning symbol indicates the level of danger.

WARNING

Indicates a dangerous situation which could lead to death or serious injury if not avoided.

CAUTION

Indicates a dangerous situation which could lead to minor injuries if not avoided.

Useful tips, recommendations and information to ensure efficient and smooth operation.

This product complies with the directive 93/42/EEC for medical products. The launch date for this product is specified in the CE declaration of conformity.

Manufacturer

1.3 Intended use

This product is intended to be used solely as a commode. Any other use can cause injury or shorten product service life.

WARNING!

Risk of injury

– Do not use this product in the shower.

1.4 Scope of delivery (Fig. 1)

Ⓐ	Frame	Ⓕ	Seat
Ⓑ	Armrest tube	Ⓖ	Pan
Ⓒ	Armrest pad	Ⓗ	Leg
Ⓓ	Backtube	Ⓖ	Foot
Ⓔ	Lid		

1.5 Warranty

We provide a manufacturer's warranty for the product in accordance with our General Terms and Conditions of Business in the respective countries. Warranty claims can only be made through the provider from whom the appliance was obtained.

1.6 Standards and regulations

Quality is of the utmost importance for our company; all processes are based on the ISO 9001 and ISO 13485 standards.

This product features the CE mark, in compliance with the Medical Device Directive 93/42/EEC Class I.

The product meets the requirements of the standards DIN EN 12182 (Technical aids for disabled persons - general requirements and test methods), DIN EN ISO 14971 (Medical devices - Application of risk management to medical devices) and DIN EN ISO 10993-1 (Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management system).

Invacare® is continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum. We only use REACH compliant materials and components.

For further information please contact Invacare® in your country (for addresses, see back page of this manual).

1.7 Service life

The expected service life of this product is five years when used daily in accordance with safety instructions and intended use, stated in this manual. The effective service life can vary according to frequency and intensity of use.

2 Safety

2.1 Safety information

WARNING!

– Do not use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as owner's manuals, service manuals or instruction sheets supplied with this product or optional equipment. If you are unable to understand the warnings, cautions or instructions, contact a health care professional, dealer or technical personnel before attempting to use this equipment - otherwise, injury or damage may occur.

WARNING!

– Invacare® products are specifically designed and manufactured for use in conjunction with Invacare® accessories. Accessories designed by other manufacturers have not been tested by Invacare® and are not recommended for use with Invacare® products.

WARNING!

– Snap buttons must be fully engaged and aligned in the proper height adjustment hole of each leg before use. This ensures that the legs are securely locked and an even height adjustment is achieved.

– Make certain all screws, nuts and/or bolts are tight at all times.

– Inspect rubber tips on the leg extensions for rips, wear or if they are missing. Replace them immediately if any of these conditions exist.

– Leg extensions with rubber tips must be in contact with the floor at all times.

– Toilet seat must be down before attempting to sit on the commode. Pail holders are not constructed to support the weight of an individual, only the pail and its contents. Use the arm supports only for assistance. Do not attempt to use the arms to support full body weight.

– Arms should be unlatched and dropped to assist in transfer onto/off the device. Users with limited physical or mental capabilities should be supervised or assisted when using the commode.

WARNING!

Risk of falling!

– Do not stand on this product.
– Do not use as a climbing aid.

WARNING!

Risk of injury!

– Do not use if defective.
– Do not make any unauthorized alterations or modifications to the product.

2.2 Labels and symbols on the product (Fig. 2)

Ⓐ	Conformity mark
Ⓑ	Production date
Ⓒ	Serial number
Ⓓ	Load capacity
Ⓔ	Note (observe accompanying documents)
Ⓕ	Product designation

-  The identification label is attached to the backside of the commode.

**WARNING!****Weight Limitation**

- The Commode listed has a weight limitation of 113 kg.
 – Consult a qualified technician if assistance is required in selecting the appropriate Invacare® commode for you.

3 Setup and usage

3.1 Installing the pail (Fig. 3)

1. Install the pail (A) onto the pail holder (B) in the commode frame (C).
2. Turn the pail so the tabs (D) keep the pail on the pail holder.

3.2 Installing/Removing the commode seat and lid (Fig. 4)

**CAUTION!****Risk of product damage**

- Cold seat clamps may break when removing/installed onto commode back frame rail.
 – Before removing/installing seat and lid, allow seat/ring to reach room temperature if exposed to cold.
 – Always remove the lid before removing the seat.
 – Pull up on the back corner of the lid or seat when removing.

1. Align the seat clamps (A) of the toilet seat (B) onto the back frame rail (C).
2. Push down on the seat until it clamps firmly into place.
3. Repeat steps 1-2 to install the lid (D).
4. Reverse steps 1-3 to remove the lid and seat.

3.3 Installing the backtube (Fig. 5)

1. Press in and hold both snap buttons (A) on the backtube (B).
2. While holding the snap buttons, insert both backtube endings into the frame tube supporters (C).
3. Push both backtube endings in until the snap buttons are fully engaged into both holes in the frame tube. There will be an audible click.

3.4 Adjusting the leg height (Fig. 6)

1. Align the snap button (A) of each leg frame (B) with the appropriate adjustment hole (C) in the leg extension (D).
2. Make certain that the snap buttons are fully engaged (E) into the adjustment hole of each leg extension. There will be an audible click.
3. Make sure all leg extensions are at the same height before using.

3.5 Operating the drop arms (Fig. 7)

1. Press and hold the button (A) on the armrest (B).
2. While holding the button, move the armrest all the way up or down.
3. Let go of the button and make sure armrests are locked in position.

4 Maintenance

4.1 Cleaning and lubrication

Cleaning

**CAUTION!****Risk of product damage.**

- Do not use an abrasive detergent, cleaner or cloth to clean or dry the commode.
 – Regular cleaning and inspection is recommended.

1. Clean with a damp cloth and non-abrasive cleanser.

Lubrication

To make raising or lowering the seat/lid easier, you may lubricate the commode back frame rail with a very small amount of petroleum jelly or furniture polish.

5 After Use

5.1 Re-use

The product is suitable for repeated use. The number of times it can be used depends on how often and in which way the product is used. Before reuse, clean the product thoroughly (® 4.1 Cleaning and lubrication, page 4), check the safety functions and reset if necessary.

**WARNING!****Risk of Injury**

- Always have a authorized dealer inspect the product for damage before transferring the product to a different user. If any damage is found, Do not use the product. Otherwise injury may occur.

5.2 Disposal

The disposal and recycling of used devices and packaging must comply with the applicable local legal regulations.

6 Technical data

6.1 Dimensions and weight (Fig. 8)

(A)	820 – 920 mm	(E)	380 – 480 mm
(B)	470 mm	(F)	220 mm
(C)	530 mm	(G)	360 mm
(D)	450 – 550 mm	(H)	430 mm

Product weight: 8.8 kg

Operating temperature: 10 – 40 °C

6.2 Materials

Tubes / frame / pin	Steel
Seat / pan	PP
Armrest pads	PE
Plug / Plastic rings on tubes	LDPE
Foot	PVC
Snap button	Cu/65Mn

All product components are corrosion resistant and latex free.

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Invacare Corporation
Aquatec Operations GmbH