


CE



DSL-105 proximity head control

DXTM THE
ULTIMATE
POWERCHAIR
CONTROL
SOLUTION



dynamicTM 

DSL-105

USER INSTRUCTION SHEET

End User Instructions

To keep yourself and others safe while using your DSL-105 Proximity Head Control, please observe the following cautions:



The Proximity Head Control cannot be used as a headrest or head support while the control is selected to drive the chair, because the chair will start to drive if the head comes close to the pads.

- **To use the head control as a headrest, turn off the chair first.**



The sensors of the head control are sensitive to water. This means that if enough water is present close to the sensors, they may activate and the chair may start to drive.

- **Do not operate the head control with wet hair**, with wet headwear or neckwear, in wet weather, or in any circumstance where water may come close to the sensors.



The sensor pads are made of water resistant vinyl. This material helps any water or sweat to run quickly off the pads before the sensors can activate.

- **Do not cover the sensor pads** with any material (for example, cloth for increased comfort). If the pads are covered with material that can absorb water, the water will stay on the pads, the chair will start to drive and it will not stop driving until the water has dried out.
- Regularly check the vinyl pad covering of the head control for cracks, damage or any form of wear and tear. **If the vinyl covering is damaged, do not use the head control** and have the coverings replaced immediately. Holes in the vinyl may allow water entry into the pads, which may activate the sensors and cause the chair to start driving.



The distance sensors are sensitive to electric waves, such as those that are transmitted by portable communication devices and electronic security gates.

- Do not use the DSL-105 near any device that transmits strong electric waves, otherwise the chair may start to drive unexpectedly.



If the chair drives unexpectedly, **stop the chair by pressing**

- **The on/off button**, or
- The 'Mode' button (only in 3Q mode).

Note: This instruction sheet is for guidance only and is not intended to replace the mandatory User Training.