



**Yes, you can.®**

Invacare nv  
Autobaan 22 - 8210 Loppem - België - www.invacare.be

Juillet 2021

**Invacare® Top End®  
Pro-2 All-Sport Aluminum  
Wheelchair**

CE

Date :

Référence :

Revendeur :

Veuillez indiquer votre choix en cochant les cases     Option standard     Option sans supplément

CODE	Description	Prix public €, TVA de 6% incluse
◆ TE10014	Top End Pro-2 All Sport	3250,00
<b>Sizing:</b>		
<b>Seat Width:</b>		
TE15013	14" Seat Width - Frame Color - White	
TE15017	16" Seat Width - Frame Color - Red	
TE15021	18" Seat Width - Frame Color - Blue	
<b>Seat Depth:</b>		
◆ TE15134	Adjustable Seat Depth 13" to 17"	Standard
<b>Seat Height:</b>		
◆ TE15096	Adjustable front seat height: 18" to 20" Adjustable rear seat height: 15" to 20"	Standard
<b>Back Height:</b>		
◆ TE15136	Adjustable back height 11" to 14"	Standard
<b>Wheelbase: (Front to Rear Axle)</b>		
◆ TE75068	17" Wheelbase	Standard
<b>Wheel Camber Angle:</b>		
◆ TE45008	15 Degree Camber	Standard
<b>Rear Wheel Options:</b>		
<b>Size and Style:</b>		
◆ TE50191	25"	Standard
◆ TE50196	High Performance	Standard
<b>Handrim Size/Material:</b>		
◆ TE50297	6tabs, alu	Standard
<b>Handrim Style/Mounting:</b>		
◆ TE50308	Standard Mount handrim	Standard
<b>Tires:</b>		
◆ TE50216	25" Primo, Grey tire	Standard
<b>Center of Gravity: (front of backrest to center of axle)</b>		
◆ TE45073	Adjustable 2.5" to 5" center of gravity	Standard
<b>Back Type:</b>		
◆ TE15159	Angle Adjustable, Non-Folding backrest	Standard
<b>Side Guards:</b>		
◆ TE20014	Screw-on, Aluminum sideguards	Standard
<b>Anti-Tip:</b>		
◆ TE20030	Single telescoping, removable, adjustable Height anti-tip	Standard
<b>Footrest:</b>		
◆ TE60003	Clamp-on tubular; height adjustable footrest	Standard
<b>Upholstery:</b>		
<b>Material/Color:</b>		
◆ TE30001	Nylon, Black upholstery	Standard
<b>Seat Style:</b>		
◆ TE30019	Screw-On seat upholstery	Standard
<b>Back Style:</b>		
◆ TE30012	Adjustable Tension back upholstery	Standard
<b>Options:</b>		
◆ TE30005	2" Soft Foam Seat Cushion	Standard
◆ TE25060	Hook and Loop Positioning Strap	Standard
TE30006	3" Soft Cushion	91,00



	CODE	Description	Prix public €, TVA de 6% incluse
	TE30008	2" Hard Cushion	108,00
	TE30009	3" Hard Cushion	108,00
	TE30010	4" Hard Cushion	108,00
<b>Options and Accessories:</b>			
<b>Safety:</b>			
◆	TE25066	Calf Strap	Standard
◆	TE25065	Rear Frame Protector	Standard
<b>Positioning:</b>			
	TE60030	Click strap single pad for footrest	123,27
	TE25001	Click Strap set, 8" Small	123,27
	TE25002	Click Strap set, 10" Medium	132,03
	TE25003	Click Strap set, 12" Large	138,02
<b>WINGS:</b>			
	TE 25074	Clamp-on Basketball Wing set	600,00