



INVACARE® TOP END®
FORCE™ CC
HANDCYCLE



The **Invacare® Top End® Force™ CC Handcycle** is a handcycle you can ride on wide trails that is lightweight but also includes extreme climbing gears, sturdy 559-26" bike wheels equipped with mountain bike tires and a stop-on-a-dime front disk brake. The design is a cross between the popular Invacare® Top End® Force™ RX Handcycle and the Invacare® Top End® Force™ G Handcycle with enhanced ground clearance and an aggressive back angle adjustment that puts you in charge of rough terrain. Now you can force the pace....cross country especially with the BionX® Power Boost option which adds 25%, 50%, 75% or 100% power assist.


Features

Ovalized 7000 series aluminum frame
Front disk brake
Adjustable seat length
Adjustable back angle
Aluminum fork
Adjustable crank housing
Multi-position adjustable footrests
Safety straps
Invacare® Top End® V crankset

Ovalized aluminum ergonomic handpedals
Shimano® 27 speed drive train
Shimano® Rapid Fire™ hands-on shifter brake
11-34 SRAM® cassette
Cable driven upper derailleur shifter
Upper chainrings: 44-34-22
Steering dampener
559-26" high performance wheels

Knobby or Cruiser tires
Threaded axles
Parking brake
Seat and back cushions
Chain guard
Safety flag
250 lb. weight capacity
Stealth Black powder coat paint
BionX® Power Boost option adds 25%, 50%, 75% or 100% power assist

Model #: FRCCC



1. Front disk brakes
2. BionX® Power Boost (available in the US and Canada only)
3. Optional dual rear disk brakes
4. Hands on shifting and braking
5. Welded footrest clamps
6. Brake lever for optional rear disk brakes



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

©2014

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Shimano and Rapid Fire are a registered trademark of Shimano Inc. SRAM is a registered trademark of SRAM Corp. BionX is a registered trademark of BionX International Corporation. Specifications are subject to change without notification. Form No. 12-300 Rev. 02/14 140190

Options

Heart rate monitor
Alignment gauge
Backpack hydration system
Rear disk brakes
Water bottle and cage
Cordless computer
Helmet
Safety lights, mirror
Additional safety flag
Click straps

Bike rack
Indoor training roller
Extra cable kit
Leg guard attachment
Tool and tire repair kit
Wheel bag
Handcycle travel bag
Top End Bike Jersey
Crank width adapter set
BionX® Power Boost

