


INVACARE® TOP END® FORCE™ G HANDCYCLE





The Invacare® Top End® Force™ G Handcycle design offers an adjustable, trunk-powered, erect position for low paras, amputees and competitors who can "gut" it out using abdominal, arm, shoulder and chest muscles. The Invacare Top End Force G Handcycle features an internal rib center frame design to make it super stiff. The adjustable back angle (55°-90°), back height, leg rests, crank height and fore/aft seat position allows you to dial in your ultimate power position. Are you ready to force the pace?

Features

Custom skeletal reinforced aluminum frame

Adjustable seat length Adjustable back angle

Aluminum oversized rake fork

Adjustable crank housing

Multi-position adjustable footrests

Safety straps

Invacare® Top End® V or S crankset

Ovalized aluminum ergonomic

handpedals

Shimano® 27 speed drive train

Shimano Rapid Fire[™] hands-on shifters

brake

11-34 SRAM® cassette

Cable driven upper derailleur shifter

Upper chainrings: 48-39-24

Steering dampener

26" high performance wheels

High pressure clincher tires

Threaded axles

Drafting bumper

Parking brake

Seat and back cushions

Chain guard

Safety flag

250 lb. weight capacity

Model #: FRCG





- 1. Adjustable back angle 55°-90°
- 2. Drafting bumper
- 3. Cable operated crank mounted shifter for upper chainrings
- 4. Shimano® Rapid Fire™ Shifter/brake mounted on a handpedal with optional power plate
- 5. Fork with diagonal and adjustable height and chain guard
- 6. 26" High performance wheels equipped with high pressure clincher tires and treaded axles

Options

Carbon fiber wheels

SRM® Training system

SRAM® Grip or Trigger shifters

SRAM Cassette choices:

11-23 or 11-26

53-40-28 Chainring

Continental® Grand Prix clincher tires

Continental Sprinter tubular tires

Heart rate monitor

Alignment gauge

Backpack hydration system

Disk brake

Large or small hand pedals

Power plates

Water bottle and cage

Cordless computer

Helmet

Safety lights, mirror

Additional safety flag

Click straps

Crutch holder

Bike rack

110.4

Invacare Corporation

www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

Top End

www.topendwheelchair.com Sports and Recreation Products 4501 G3rd Circle North Pinellas Park, Florida 33781-5914 (727) 522-8677 (800) 532-8677

©2014

Fax: (727) 522-1007

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols [™] and [®]. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Shimano and Rapid Fire are registered trademarks of Shimano Inc. SRM is a registered trademark of SRM Service Center Inc. USA. SRAM is a registered trademark of SRAM Corp. Continental is a registered trademark of Continental Tire Corp. Specifications are subject to change without notification. Form No. 12-450 Rev. 12/14

Indoor training roller
Extra cable kit
Leg guard attachment
Extra crank arms
Tool and tire repair kit
Wheel bag
Top End bike jersey

