



INVACARE® TOP END®

# XLT HANDCYCLE



Yes, you can.®



If recreation is your goal, the Invacare® Top End® XLT and XLT Jr. Handcycles are great fun and feature low maintenance 7 speed hubs with reverse hands-on braking. Just get out there and go! All Invacare Top End recreation handcycles are equipped with a safety flag, full composite chain guard, cushion and seat/chest restraint, as standard, for extra value, comfort and safety.

---

The **Invacare® Top End® XLT Handcycle** is designed with the recreational handcycling enthusiast in mind. Equipped with seven speeds, the XLT handcycle will easily handle most terrain, whether it's riding around town or taking a spin in the country. The lightweight, oval tube aluminum XLT frame with chrome-moly steel fork features simple geometry for maximum stiffness and easy cornering. The seat and back are totally adjustable for a comfortable fit. The narrow rear-end enhances aerodynamics and the 15° of camber increases stability. If you want a great way to exercise, cross-train or just have fun, choose the XLT handcycle.

---

**Model #: XLT**



1. Twist Shifter
2. Quad Twist shifting



---

**Features**

Choice of 15" wide or 17" wide seat widths  
6061 Aluminum frame/chrome-moly fork  
Adjustable seat, back angle and crank height  
Full chain guard, chain tension idler and safety flag  
Adjustable footrest system with safety straps  
Seat and back cushions with chest restraint  
Shimano® 7 speed hub with reversing drum brake  
175 mm crankset with one chain ring  
Quad Twist shifter option helps those with hand impairments

Frictionless ovalized ergonomic hand pedals  
Parking brake  
26" High Performance wheels  
High pressure clincher tires  
Quick release axles  
15° camber  
Road crown compensator  
Reflectors  
250 lb. weight capacity  
Many options and aftermarket accessories are available



---

The **Invacare® Top End® XLT Jr. Handcycle** is a recumbent handcycle built just for kids. Easy to operate, this handcycle has seven speeds and an aluminum frame with a chrome-moly steel fork. The enhanced aerodynamics and the 15° of camber increase stability. A tall, narrow back offers maximum range of motion and support, with a chest strap for safety. This extremely stable and maneuverable three-wheel handcycle is great fun for anyone who's between 38" to 60" tall. Standard safety features include: reflectors, full chain guard, footrest safety straps, safety flag, cushion and seat restraint. The XLT Jr. Handcycle is totally adjustable and ready to ride!

---

**Model #: XLTJR**



1. Parking brake
2. Optional horizontal handpedals
3. Ergonomic handpedals

---

**Features**

6061 Aluminum frame/chrome-moly fork  
14" wide adjustable seat, back angle and crank height  
Full chain guard  
Chain tension idler and safety flag  
Adjustable footrest system with safety straps  
Seat and back cushions with chest restraint  
Shimano® 7 speed hub with reversing drum brake  
155 mm crankset with one chain ring  
Vertical ergonomic handpedals

Parking brake  
20" cruiser wheels and tires  
Quick release axles  
15° camber  
Road crown compensator  
Reflectors  
250 lb. weight capacity  
Many optional accessories are available

**Invacare Corporation**

[www.invacare.com](http://www.invacare.com)

**USA**

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Canada**

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

**Top End**

[www.topendwheelchair.com](http://www.topendwheelchair.com)  
Sports and Recreation Products  
4501 63rd Circle North  
Pinellas Park, Florida  
33781-5914  
(727) 522-8677  
(800) 532-8677  
Fax: (727) 522-1007

©2014

Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and ©.  
All trademarks are owned  
by or licensed to Invacare  
Corporation unless  
otherwise noted. Shimano  
and Rapid Fire are registered  
trademarks of Shimano Inc.  
Form No. 12-510  
Rev. 02/14  
140187

