



Yes, you can.®

Invacare nv  
Autobaan 22 - 8210 Loppem - België - www.invacare.be

Invacare® **Top End®**  
**Li'l Excelerator Handcycle**



Date :

Référence :

Revendeur :

Veuillez indiquer votre choix en cochant les cases     Option standard     Option sans supplément

**CHAIR TYPE**

LXC    Top End Li'l Excelerator Handcycle    3600.00

Seat Width: 14"

Frame Color: Royal Blue

**OTHER THAN STANDARD FRAME COLOR:**

P3    159.42

**CONTEMPORARY COLORS**

Wet Black    White Snow    Silver Metallic  
Stealth Black

**WAY COOL COLORS**

Midnight Blue    Grasshopper Green    Electric Blueberry


**TOO HOT COLORS**

Tangerine Red    Raspberry Red    Grape Madness  
Safety Yellow    Orange    Bubblegum Pink

**SAFETY AND COMFORT OPTIONS:**

	Seat Cushion	Standard
	Chain Guard	Standard
	Safety Flag	Standard
	Seat restraint: hook and loop	Standard
XLT139	Quad twist shifter option	62.60
XCL 3	Horizontal foam covered handpedals	42.38
XLT12	Mirror	11.59
XLT13	Water Bottle & Cage	40.96
XLT15	Helmet: choose size	58.12
	Small/ Medium    Medium/Large	
XLT6	Front Safety Light	20.87
XLT7	Rear Flashing Light	20.87
XLT17	Crutch Holder & Strap	105.13
XLT20*	Bike Rack for Trailer Hitch	546.63
	1 1/4" (Standard)    2" (No Charge)	
XCL13	Seat restraint: Auto style fastener	53.61

\* Indicates an item where the standard discount does not apply or there are additional shipping charges or both. Please request a quote.



Includes:

- Carbon steel frame and fork
- Adjustable seat: forwards/backwards
- Adjustable back height/ Fold over back upholstery
- Adjustable crank height
- Chain guard
- Chain tension idler
- Adjustable height/angle footrest system
- Shimano 7 speed hub w/reversing drum brake
- 155 mm crankset w/one chain ring
- Oval Ergonomic handpedals
- Parking brake
- 20" Cruiser wheels and tires
- Quick release axles
- 9° camber
- Reflectors
- Cushion, seat restraint (hook and loop) and safety flag

Frame – 1 3/4" x .065 ovalized carbon steel tubing.

Frame Design – Single main tube, upright position.

Wheelbase – Overall length 56"

Weight – Overall weight 45 lbs.

Turning Radius – 10 feet for 360° circle, 7 feet w/one backing maneuver.

Static Tip Angle - 30° for 150 lb subject

Transfer Seat Height – (w/2" Cushion) 17 1/2".

Crankset – Sugino 155mm w/chainring, adjustable height. Also includes chain tension idler.

Chainguard- Full length shock resistant composite construction.

Standard Propulsion System - -155mm Crankset w/one chainring, adjustable height crank, Shimano 7 speed front hub w/internal hub brake and 7 speed twist shifter.

Braking System – Shimano internal drum brake (coaster brake) and Chain, activated by pedaling backwards.

Gear Shifter- Shimano Tourney revoshifter indexed twist grip with gear display. For customers with limited hand function the Quad twist shifter option, a 3 pronged adaptation that is secured around the twist shifter is recommended.

Parking Brake – Caliper brake w/locking brake lever.

Handpedals – Oval Ergonomic handpedals. Optional horizontal style is available.

Options for impaired hand function–The quad twist shifter adaptation option will allow the users to easily change gears.

Hub – Front: Shimano 7 speed hub w/internal hub brake, freewheel back up feature. Rear: precision, heat treated, black, anodized w/ 1/2" quick release, stainless axles.

Performance Wheels – Rear hubs are precision heat treated, silver anodized w/quick release stainless axles. Quick release axles are easiest for transporting the handcycle. 20" x 1.75" treaded Tire,

20 x 1.5" Tube 40 psi,

20" Spoke Wheels,

20" Spoked Front Wheels

Camber - 9° Camber for stability.

Gears –7 Speeds

Seating/Back Height – 14" width, multi position for and aft sliding seat system. Back Height adjustable, angle is 100°.

Upholstery – Nylon w/hoop and loop adjustable tension seat sling and fold-over seat back. Hook and loop seat restraint and foam seat cushion w/cover included.

Footrest – Adjustable height and angle, dual footrests w/safety straps.

Warranty – 3 year on frame to original owner /limited parts warranty.

Assembly – Shipped completely assembled except for rear wheels. Shipping weight is 70 lbs.

Options:

Helmet

Quad twist shift adaptation

Seat Restraint: auto style

Water Bottle and Cage

Front and Rear Safety Lights

Crutch Holder and Strap

Mirror

Bike Rack