

# Invacare® Top End® Force CC™ Handcycle

## Ireland Dealer Order Form January 2014

For accuracy of ordering, please fax order to Invacare Customer Service on 003531 8107085 or e-mail to ordersireland@invacare.com



www.topendwheelchair.com

Dealer Name _____	Date _____	Top End Job Number _____
Account No _____	PO No _____	Phone No _____ Fax No _____
Address _____	Post Code _____	
Customer Name _____	Phone No _____	
Address _____	Post Code _____	
<b>Ship order to:</b> <input type="checkbox"/> <b>Provider</b> <input type="checkbox"/> <b>Customer</b>		

**Order will ship to Provider unless otherwise noted.**

### TOP END FORCE CC HANDCYCLE:

**FRCCC Top End Force CC Off Road 27 Speed Handcycle** €4,995

#### SEAT WIDTH: (Select One)

- 15" Stealth Black Frame and Fork**
- 17" Stealth Black Frame and Fork**

#### BRAKING SYSTEM: Front Disk Brake and Parking Brake..... **Standard**

**XLT88 Dual Rear Disk Brake upgrade.....**€325

#### CRANKARMS: 170mm length x 14" wide..... **Standard**

170mm length x 17" wide..... **No Charge**

#### WHEELS: 559-26" wheels with threaded Axles


#### TYRES:

- Knobby mountain bike tyres..... Standard**
- Cruiser clincher tyres..... No Charge**

### SAFETY & COMFORT OPTIONS

**NOTE: All non-standard options will be shipped separately**

- |                                              |                                                                                  |                                                  |
|----------------------------------------------|----------------------------------------------------------------------------------|--------------------------------------------------|
|                                              | Seat and Back Cushions .....                                                     | <b>Standard</b>                                  |
|                                              | Chain Guard.....                                                                 | <b>Standard</b>                                  |
|                                              | Safety Flag .....                                                                | <b>Standard</b>                                  |
| <input type="checkbox"/> <b>XLT9</b>         | Leg Guard Attachment .....                                                       | €62                                              |
| <input type="checkbox"/> <b>XLT84</b>        | Crank width adaptor set (adds 1.5" per set) .....                                | €20                                              |
| <input type="checkbox"/> <b>XLT71</b>        | Tool and Tyre Repair Kit .....                                                   | €70                                              |
| <input type="checkbox"/> <b>XLT12</b>        | Mirror, Leg Rest Mount .....                                                     | €24                                              |
| <input type="checkbox"/> <b>XLT13</b>        | Water Bottle & Cage .....                                                        | €29                                              |
| <input type="checkbox"/> <b>XLT28</b>        | Backpack Hydration System .....                                                  | €83                                              |
| <input type="checkbox"/> <b>XLT141</b>       | Heart Rate Monitor .....                                                         | €83                                              |
| <input type="checkbox"/> <b>XLT14</b>        | Cordless computer .....                                                          | €62                                              |
| <input type="checkbox"/> <b>XLT15</b>        | Helmet .....                                                                     | €62                                              |
| <input type="checkbox"/> Small (6 5/8" – 7") | <input type="checkbox"/> Medium (6 7/8" – 7 1/2")                                | <input type="checkbox"/> Large (7 1/4" – 7 5/8") |
| <input type="checkbox"/> <b>XLT16</b>        | Extra Safety Flag .....                                                          | €24                                              |
| <input type="checkbox"/> <b>XLT6</b>         | Front Safety Light .....                                                         | €20                                              |
| <input type="checkbox"/> <b>XLT7</b>         | Rear Flashing Light .....                                                        | €20                                              |
| <input type="checkbox"/> <b>XLT 33</b>       | Alignment Gauge .....                                                            | €62                                              |
| <input type="checkbox"/> <b>XLT40</b>        | Cable Kit.....                                                                   | €31                                              |
| <input type="checkbox"/> <b>XLT74</b>        | Extra Steering Dampener .....                                                    | €13                                              |
| <input type="checkbox"/> <b>XLT42</b>        | Indoor Stationary Trainer .....                                                  | €214                                             |
| <input type="checkbox"/> <b>TER126</b>       | Click Strap Large (pad length~12") .....                                         | €72                                              |
| <input type="checkbox"/> <b>TER125</b>       | Click Strap Med (pad length~10") .....                                           | €67                                              |
| <input type="checkbox"/> <b>TER124</b>       | Click Strap Small (pad length~8") .....                                          | €63                                              |
| <input type="checkbox"/> <b>RACE22</b>       | Top End T-Shirt <input type="checkbox"/> L <input type="checkbox"/> XL .....     | €10                                              |
| <input type="checkbox"/> <b>RACE4</b>        | Wheel Bag .....                                                                  | €315                                             |
| <input type="checkbox"/> <b>RACE20</b>       | Handcycle Travel Bag .....                                                       | €350                                             |
| <input type="checkbox"/> <b>RACE13</b>       | Top End Bike Jersey <input type="checkbox"/> L <input type="checkbox"/> XL ..... | €78                                              |



### FORCE CC HANDCYCLE SPECIFICATIONS

- Frame** – Ovalized 7000 series aluminium, Stealth Black frame colour
- Wheel base:** 60" Standard, Overall Length 86"
- Removable Fork** – 7000 series aluminium with adjustable crank housing.
- Components:** Shimano® Deore Rapid-fire® Brake/Shifter
- Mounting Position: Shifter / Brake Lever:** Right Handle
- Front Wheel Cassette:** 11-34
- Chain ring:** 44-34-22
- Handpedal Size:** Medium (1 3/4")

#### Diagonal and Vertical Crank Height Adjustment

**Camber:** 0 degrees

**Turning Radius** – 18 feet for 360° circle, 12 feet with one backing manoeuvre.

**Static Tip Angle** - 39° for 250 lb/113kg subject.

**Transfer Seat Height** – 11.5" with 2" cushion.

**Ground clearance-** 6.5"

**Backrest** – Adjustable fore and aft.

**Back Angle** – Adjustable angle back of 90° - 120°.

**Propulsion System** – Shimano components 27 speeds, 44-34-22 chain rings and 11/34 cassette. The 11/34 cassette is operated by the Hands-On Rapid Fire unit which is mounted on the right hand pedal. It operates both the 11/34 cassette and the hands-on front disk brake. The top chain rings utilizes a cable driven shifter system located on upper crank.

*All prices subject to change without notice.  
Subject to Invacare's standard terms and conditions.*

## Invacare® Top End® Force CC™ Handcycle

**Braking System**-A parking brake/second brake is mounted on the upper fork. Front disk brake is standard and the brake lever is mounted on the right handpedal. Rear disk brakes are optional and the brake lever is mounted on the left side of the fork.

**Crank set** – Top End Crankset

**Steering** – Pivot steering and raked fork for controlled turning and stability.

**Crank** – Adjustable height both vertical and diagonal.

**Hand pedals** – Hand pedals are made of lightweight aluminium and are ovalized to precisely fit the hand. The bearings allow the pedal to rotate while pedaling and are located in the crank set and not in the pedal, which makes this pedal very comfortable.

**Footrest** – Multi-adjustable fore and aft with straps to securely cradle calves and feet. Straps are 4" wide and feature safety straps for ankles.

**Performance wheels** –Silver spoked wheels with threaded axles and alloy rims. 14 gauge stainless steel spokes. 50-559 (26)" x 1.95" Knobby Clincher Tyre/Tube 40- 60psi. or optional 26" x 1.5"Cruiser Clincher Tyre.

**Upholstery** – Black nylon upholstery with 2" seat and back cushion standard. Seat restraint included.

**Warranty** –3 year on frame and fork to original owner/limited parts warranty.

**Weight Limit:** 250lbs/113kg

**Shipping weight:** Shipping weight is 62lbs/28kg Box size is 77" x 27" x 23".



**Invacare Ireland Ltd** [www.invacare.ie](http://www.invacare.ie)

Unit 5 Seatown  
Business Campus  
Seatown Road  
Swords

T: 003531 8107084  
F: 003531 8107085

---

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Shimano and Rapid Fire are registered trademarks of Shimano Inc. LeMond and Revolution are a registered trademarks of Lemond Fitness, Inc.  
©2014 Invacare Corporation