



INVACARE® TOP END®  
**FORCE™ -3**  
HANDCYCLE



The **Invacare® Top End® Force™ -3 Handcycle** is a great handcycle at a great price. This bike is ready to ride with 27 speeds, performance wheels and Shimano® components. Perfect whether you are just looking to get started in handcycling or ready to force the pace at your next competition.

---

## Features

15" or 17" seat width

Stealth black powder coat paint

60" Wheelbase

Adjustable back angle, seat and footrests

27 speed drive train

Shimano® Deore Rapid Fire™ shifter/brake

Shimano® 170 mm cranks

Avid® disk brakes

Shimano® cable driven upper derailleur

Adjustable crank height: diagonal and vertical

Alex DM18 26-559 wheels

0° camber

High pressure cruiser clincher tires

Seat and back cushions

Chest and footrest safety straps

Drafting bumper

Threaded axles

Cable springs

Frictionless ovalized ergonomic hand pedals

Safety flag

Chain guard


Reflectors

Tektro® parking brake

250 lb. weight capacity

---

## Model #: FRC3



1. Wraparound frame with adjustable seat upholstery
2. Cable driven shifting for upper chainrings
3. Front disk brake and welded footrest clamps
4. Backrest adjusts forward and backward
5. Adjustable back angle allows 35-65° recline
6. Shimano components for hands-on shifting and braking

---

## Options

Heart rate monitor

Alignment gauge

Backpack hydration system

Water bottle and cage

Cordless computer

Helmet

Safety lights, mirror

V-Crank upgrade

Trunk power position backrest

BionX® Power Boost

Additional safety flag

Click straps

Bike rack

Indoor training roller

Knobby tires

Extra cable kit

Tool and tire repair kit

Wheel bag

Top End Bike Jersey

Factory assembly

---

## Packages

### **Invacare® Top End® Force™-3 Cross Country Off Road Handcycle Package**

can be ordered with a trunk-powered position short back, aggressive knobby tires and a V crankset that has a 44 – 34 – 22 chain ring to help you power up the hills. The 26 – 559 wheels can accept both cruiser or knobby tires! This package will take the place of the Invacare® Top End® Force™ CC Handcycle and offers great value.

### **Invacare® Top End® Force™-3**

**Trunk Power Handcycle Package** is configured by adding the straight back style/ trunk- powered position back, the standard 26 – 559 cruiser tires and the V crankset with chain rings that are 48 – 39 – 24.

### **Invacare® Top End® Force™-3 with BionX® Power Boost Handcycle Package**

is a great option for people who need just a bit of help for climbing hills or getting back from a long ride. The BionX power boost can provide 25%, 50%, 75% or 100% assist.



**Invacare Corporation**

[www.invacare.com](http://www.invacare.com)

**USA**

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Canada**

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

**Top End**

[www.topendwheelchair.com](http://www.topendwheelchair.com)  
Sports and Recreation  
Products  
4501 63rd Circle North  
Pinellas Park, Florida  
33781-5914  
(727) 522-8677  
(800) 532-8677  
Fax: (727) 522-1007

©2014

Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and ®.  
All trademarks are owned  
by or licensed to Invacare  
Corporation unless  
otherwise noted. Shimano,  
Rapid Fire and Tektro are  
registered trademarks  
of Shimano Inc. Avid is  
a registered trademark  
of SRAM Corporation.  
BionX is a registered  
trademark of BionX  
International Corporation.  
Specifications are subject  
to change without  
notification.  
Form No. 13-692  
Rev. 12/14  
141215



**Yes, you can.®**