

Orion[®] Metro

Feel the Freedom



Quality, comfort and safety

for independence, confidence and style



Yes, you can.[®]

Feel the difference

Our scooters have been designed together with healthcare experts, to provide unrivalled standards of quality, safety and comfort.

Through continuous dedicated research and an unrivalled experience in designing for mobility, at Invacare we pride ourselves in scooters that improve quality of life. Being able to regain confidence, trust and independence from your scooter is important.

From ergonomic steering controls to the adaptable cushioned seats and adjustable seat risers, every detail of our scooters is considered and refined to make getting about feel great.



Quality



Market-leading build quality with exceptional attention to detail.

Safety



Compliant with the latest European Medical Device Standards.

Comfort



Adjustable design to provide lasting comfort.

Feel in control with Orion® Metro

Introducing our Orion® Metro, a high-quality scooter designed to enhance safety, comfort and added driving performance. The Orion® Metro offers everything you need from a mobility scooter.

Compact and easy to drive, this scooter handles a wide variety of surfaces. Designed with a powerful motor, the Orion® Metro drives across uneven ground, around corners and over obstacles without

compromising safety or comfort. Incredibly simple to steer and start, the Orion® Metro puts you in total control.

up to
12 km/h
max speed

up to
54 km
max distance
on a single
charge

136 kg
max user
weight



Improved driving control

- Ergonomic steering for added control
- Designed for all-terrain performance
- Anti-splash guards



Safety in mind

- Automatic speed reduction for safe cornering
- Two-step disengaging lever to prevent free-wheeling
- Brake lights and indicators
- Headlights for improved visibility



Built around comfort

- High quality cushioned seat
- Adjustable steering column
- Adjustable seat height
- Flip-up armrests and swivel seat makes getting on and off easy



Adjustable steering for greater comfort



Designed for
**improved traction
around corners**
and safe, all-terrain
performance



Optional connectivity

Choose Zeta™ connectivity, an innovative technology allowing access to your mobility scooter's battery and journey data via an easy-to-use mobile app.



Swivel seat for easy transfer



High quality lighting



Two-step disengaging lever to prevent free-wheeling



Add to your scooter

Tailor the Orion® Metro to your taste and needs. Turn over to explore the optional extras available.

Add to your scooter

We have a wide range of optional extras available, so you can customise your Orion® Metro scooter. Ask your retailer to discover the full range.



Rear basket



Side & rear view mirrors



Front & rear bumper



Mobile phone holder



Crutch holder



Seatbelt



Cupholder



Interchangeable shroud panels

Technical data

For more information about this product, including the user manual, please visit our website www.invacare.eu.com

	Backrest width ▼	Seat depth ▼	Seat to floor height ▼	Backrest height ▼	Driving unit width ▼	Total length ▼
3 wheels	510 mm	470 mm	Min. 440 mm Max. 510 mm	490 mm	660 mm	1240 mm
4 wheels	510 mm	470 mm	Min. 440 mm Max. 510 mm	Comfort: 500 mm Deluxe: 490 mm	660 mm	1270 mm
	Empty weight ▼	The heaviest part ▼	Max. user weight ▼	Battery capacity ▼	Turning diameter ▼	Max. climbable obstacle height ▼
3 wheels	103 kg	52 kg	136 kg	2 x 12 V/50 Ah (C20) AGM	2250 mm	80 mm
4 wheels	110 kg	57 kg	136 kg	2 x 12 V/50 Ah (C20) AGM	2600 mm	80 mm
	Ground clearance ▼	Max. safe slope ▼	Motor capacity ▼	Max. driving distance ▼	Max. driving speed ▼	
3 wheels	80 mm	10 °	12 km/h: 550 W / 1500 W ¹	54 km	12 km/h	
4 wheels	80 mm	10 °	6 km/h: 240 W / 500 W ¹ 10 km/h: 240 W / 600 W ¹ 12 km/h: 550 W / 1500 W ¹	54 km	6 / 10 / 12 km/h	

¹ Peak performance

Colours*



Marble White



Alium Blue



Diamond Silver



Garnet Red



Black Vinyl

* Scooter colours may vary slightly. Please check with your retailer.

Upholstery

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland

