

# Invacare Colibri®

Feel the Freedom



Stylish, lightweight and transportable

Designed with Invacare Litelock™ technology



Yes, you can.®

# Feel the difference

Our scooters have been designed together with healthcare experts, to provide unrivalled standards of quality, safety and comfort.

Through continuous dedicated research and an unrivalled experience in designing for mobility, at Invacare we pride ourselves in scooters that improve quality of life. Being able to regain confidence, trust and independence from your scooter is important.

From folding transportable designs, cushioned seat and adjustable seat risers, every detail of our scooters is considered and refined to make getting about feel great.

## Quality



Market-leading build quality with exceptional attention to detail.

## Safety



Compliant with the latest European Medical Device Standards.

## Comfort



Adaptive and adjustable design to provide lasting comfort.



# Feel in control with Invacare Colibri®

Introducing the Invacare Colibri®, our lightest and most transportable scooter. Designed to complement your style with a choice of colours and improve your independence.

Compact and easy to drive, this scooter has been designed with a patented Litelock™ system making it easy to transport into the back of your car. The swivel seat makes getting on and off your scooter simple.

The Colibri® is incredibly simple to start and steer. It has been designed to give you a comfortable and safe drive with the ability to manoeuvre in tight spaces.

up to  
**8km/h**  
max speed

up to  
**16km**  
max distance  
on a single  
charge

**136 kg**  
max user  
weight



Easy transportation with the Invacare Litelock™ system



## Easy to use and transport

- Travel easily in style
- Ergonomic steering for added control
- Simple to steer and drive



## Safety in mind

- Two-step disengaging lever to prevent free-wheeling
- Brake lights, headlights and indicators
- Additional optional extras available



## Built around comfort

- High quality cushioned seat
- Adjustable seat height
- Flip up armrests and swivel seat makes getting on and off easy



## Adjustable steering for greater comfort



Our **lightest** and most **transportable** scooter



**Easy battery charge**  
in situ or remove from scooter



Swivel seat for easy transfer



High quality lighting



Adjustable steering for improved seating position




## Add to your scooter

Tailor the Invacare Colibri® to your taste and needs. Turn over to explore the optional extras available.

## Add to your scooter

We have a wide range of optional extras available, so you can customize your Invacare Colibri® scooter. Ask your retailer to discover the full range.



Additional battery pack

Large mirror

Crutch holder

Optional size wheelsets

'Easy fit' suspension kit

Interchangeable shroud panels

## Technical data

For more information about this product, including the user manual, please visit our website [www.invacare.eu.com](http://www.invacare.eu.com)



Backrest width ▼

465 mm



Seat depth ▼

400 mm



Seat to floor height ▼

370 - 420 mm



Backrest height ▼

385 mm



Driving unit width ▼

505 mm



Total length ▼

1010 mm



Empty weight ▼

44 kg (12 Ah)  
49 kg (18 Ah)



Max. user weight ▼

136 kg



Battery capacity ▼

2 x 12 Ah  
2 x 18 Ah



Turning diameter ▼

2200 mm



Max. climbable obstacle height ▼

45 mm



Max. safe slope ▼

6° / 10.5%



Motor capacity ▼

200 W



Max. driving distance ▼

Up to  
16 km



Max. driving speed ▼

8 km/h

## Colours\*



Marble  
White



Sapphire  
Blue



Diamond  
Silver



Jasper  
Red



Jade  
Green

## Upholstery



Black Vinyl

\* Scooter colours may vary slightly.  
Please check with your retailer.

**Tel: +41 61 487 70 70**  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
Neuhofweg 51  
4147 Aesch BL  
Switzerland

