



Rea *Azalea Minor* – designed for petite adults and teenagers

The **Rea *Azalea Minor*** is designed for users that require a seat that has a smaller depth and width. Adapted from the **Rea *Azalea***, the ***Azalea Minor*** boasts all the advantages of a reliable, tilt-able wheelchair with its special weight-shifting mechanism. This model is ideal for all users that require a solid, stable wheelchair.





Features and Options

Rea Azalea compatibility continues

A new addition to the **Rea Azalea** Family, **Azalea Minor** is designed for users that need a smaller-sized seat, such as petite adults or teenagers. One of the key features of **Azalea Minor** is the 'Tilt in Space' weight-shifting mechanism, which ensures a stable tilt and makes maneuvers easier. As with all the **Azalea** models, the width and depth of the **Azalea Minor** can be adjusted to suit individual user needs.

Push handles outward angled

The push handles are angled outwards, increasing comfort for caregivers.



X-Harness

A soft harness supports the trunk and is always used in conjunction with the pelvic belt.

For other **Rea Azalea** accessories, please consult the **Rea Azalea** brochure



Laguna back cushion

The inner dimensions are as large as the **Azalea**, at 39 cm wide.



Shorter Chair

The chassis is 3 cm shorter, the seat depth is 5 cm shorter and the backrest 10 cm lower than on standard models.

Technical data

Rea Azalea Minor	34 - 44 cm	38 - 45 cm	40/45 cm	55 - 65 cm	24 - 36 cm	33 - 50 cm	-1° - 20° (22/24") / -1° - 25° (16")
Rea Azalea Minor	0° - 30°	SW + 22 cm (SW + 23 cm)	90 - 125 cm	90 - 97 cm	28 kg (SW 34) / 34 kg (SW 44 cm)		
Rea Azalea Minor	Max 75 kg	Transport weight 15 kg	Crash-tested ISO 7176-19				

Rea Azalea Minor complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information regarding vehicular transportation and for more pre-sale information for this wheelchair, please consult the following pages on our website: www.invacare-rea.com under Rea wheelchairs/Technical downloads.

Invacare International Sarl
Route de Cité-Ouest 2 - 1196 Gland Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeuropa@invacare.com
www.invacare.eu.com

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Rea Azalea Minor - EU - 08/2010

Frame colour

Upholstery



Pearl Grey



Grey Dartex TR23



Grey plush TR 32



Yes, you can.®

