

# DOLOMITE | JAZZ



## New generation...

**DOLOMITE JAZZ** is the new generation rollator that goes all the way. A stylish design with smart functions. Carefully thought-out, well designed components make **DOLOMITE JAZZ** a winner. Continuous vertical adjustment and ergonomic handle with built-in brake system that surpasses its competitors. Positioning marks on handlebars to facilitate height adjustment.





**Foldable**

Spacious basket made of durable material that is flexible when folding. Thanks to welded aluminium frames, the **DOLOMITE JAZZ** is stable and steady throughout its service life. Simple and balanced folding with the help of a seat strap. Easy to handle, occupies a small space and stands in folded position. Nice soft seat.

**Adaptable to individual needs**

**DOLOMITE JAZZ** is available in two different heights. Choose the higher seat height when having difficulties to change from sitting to standing. Easily adjustable handle that can be set at precisely the right height.

**NEW!** Slow-down brake - A safety feature that prevents the walker from rolling a way.

**Features and Options**



**Comfortable backrest**

A comfortable and safe support to lean against. Mounted in the handle bars.



**Cane holder**

Easy to fit, easy to use. No tools.



**Basket cover**

To protect what you have in the basket from rain.



**Slow-down brake**

A safety feature that prevents the walker rolling away. Perfect for gait training and one handed solutions.

**Technical data**

JAZZ 600	62 cm	67 cm	44 cm	75 - 97 cm	60 cm
JAZZ 520	62 cm	67 cm	44 cm	67 - 81 cm	52 cm
JAZZ 600	39 cm	20 cm	7,8 kg	130 kg	26 cm
JAZZ 520	39 cm	20 cm	7,2 kg	130 kg	26 cm

**Colour**



Dark gray

**Article code**

**JAZZ 600** 12140-46-24

**JAZZ 520** 12142-46-24

Including Tray and Basket

INVACARE® DOLOMITE AB  
Växjövägen 303 75  
Box 200, S-343 75 Diö  
Tel. +46 (0)476 535 00  
Fax +46 (0)476 535 99  
www.dolomite.biz

© 2010 Invacare® International Sarl  
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Jazz - EU - 09/2010



**Yes, you can.®**

