

# DOLOMITE | *MELODY*



## **Turns life into a dance**

5 lines positioning statement to include the following facts;  
The new **DOLOMITE MELODY** is a neat stylish walker with features that far outstrip the competitors in the same price class. Easy to deal with as it can be folded and stands in the folded position. The wireless braking system is hidden inside the frame. **DOLOMITE MELODY** is so neat and easy to handle it makes you wanna dance.





### Wireless, integrated brakes

No brake cable to get caught up. Components in stainless steel for long life. Patented braking mechanism with tried and tested function.

### Stands in folded position

Standing when folded makes the walker easier to park and store on a small surface. Push the seat to get the rollator to operating position.

### Easier to walk and reach things

The seat is secure and comfortable. It can be folded up to gain more space for walking and to make it easier opening doors or reaching things for example kitchen tops.

### Infinite adjustment

Using small steps, precise height and angle can be set. In this way, the walker can be customized for different arm lengths and hand positions. Combined the unique adjustment possibilities with

three height models and a wide range of accessories **DOLOMITE MELODY** is one of the most flexible rollators on the market.

## Features and Options



### Practical basket

that even holds small items.  
Max load: 10kg.



### Tray

Made of impact-resistant plastic.  
Four stable feet.



### Back strap

A comfortable and safe support to lean against. Mounted into handle bars.

## Features and Options



### Cane holder

Makes the cane easily accessible.



### Slow down brake

Ideal for gait training or when braking is a challenge.



### Kerb climber

"Pedal" that simplifies mounting kerbs or similar.



### Oxygene tank holder

Secure mounting of oxygene tube on the walker.



### One-hand brake

Both brakes are activated with one hand.  
Order complete product only.



### HP-handle set

For stroke patients.  
Order complete product only.



### Adapt handles set

Anatomic grips that can be adjusted to any position.  
Order complete.



### Drip feed holder

Holds drip feed on place. Even for outdoor use.



### Wheel guard

Protects clothing from being dirty from wheel splash.

## Technical data



Melody S	55 cm	63 cm	47 cm	66 - 81 cm	48 cm
Melody M	61 cm	65 cm	54 cm	70 - 90 cm	53 cm
Melody L	61 cm	69 cm	54 cm	78 - 98 cm	61 cm



Melody S	36 cm	20 cm	7 kg	135 kg
Melody M	43 cm	20 cm	7,2 kg	135 kg
Melody L	43 cm	20 cm	7,4 kg	135 kg

## Colour



RAL 1035

## Article code

<b>Melody S</b>	<b>12014</b>
<b>Melody M</b>	<b>12012</b>
<b>Melody L</b>	<b>12010</b>

INVACARE® DOLOMITE AB  
 Växjövägen 303 75  
 Box 200, S-343 75 Diö  
 Tel. +46 (0)476 535 00  
 Fax +46 (0)476 535 99  
 www.dolomite.biz

© 2010 Invacare® International Sarl  
 All rights reserved.

All information quoted is believed to be correct at time of print.  
 Invacare® reserves the right to alter product specifications without prior consultation.

Melody - EU - 09/2010



**Yes, you can.®**