

Accent[®]

For everyday care



Medical profiling bed

Trusted by healthcare professionals



Affordable quality

The Accent bed provides a simple and affordable quality solution for both homecare and long-term care settings.



Safety

All Invacare beds meet IEC 60601-2-52 safety standards.



Quality

Our bed designs and specifications adhere to the highest standards.



Care-led

Each Invacare bed is designed with the knowledge of care settings in mind.



Comfort

Every Invacare bed promotes rest, dignity and independence.




Homecare ready

Accent is an easy and user-friendly choice for patients, families and caregivers.



Shown with optional side rails. Invacare mattress sold separately.

Designed with care



With straightforward assembly, smooth adjustments and thoughtful details, Accent supports quality care every day.

01



Affordable choice
A simple and affordable quality bed.

02



Quick & easy setup
Simple tool-free assembly and disassembly. Easy to handle.

03



Adjustable positioning
Backrest and thigh sections adjust at the touch of a button.

04



Height adjustable
Adjusts between 40 - 80 cm for comfortable working heights and optimal patient care.

05



Intuitive hand control
Clear visual display for ease of use. Optional lockable tilt available.

06



Range of side rails
Choose from the Slim, Scala 2, Verso HC or Aria options.



Easy to clean
Suitable to wipe with mild detergents making it ideal for single or multiple use environments. See User Manual.



Compact power unit & cabling
Smaller unit with neater and detachable cables gives more space under the bed and makes moving the bed easier.

08

07

09

11° Supportive tilt
For positioning and pressure management.



Accent® options

Simplifying care with a choice of everyday options.



Lifting pole



Slim
Full length side rail in Pine wood.



Scala 2
Side rail in steel.



Verso HC
Collapsible side rail in steel.



Aria
Full length side rail in steel.



Bed end
Optional panel.



Alternative handset
With lockable tilt for peace of mind.

Easy to transport

Optional Transportation Kit makes moving and storing the bed even easier.



Watch the How to video on assembly & disassembly.



Tech Support & After Sales

Ask our team. We are always here to help.

Tel: +41 61 487 70 70
support@invacare.com

Technical specifications



Width inside
90 cm



Width outside
102 cm



Height adjustment
40 - 80 cm



Length inside
208 cm



Length outside
224 cm



Tilt angle
11°



Max user weight
145 kg



Total weight
64 kg



Heaviest part
16 kg



Mattress support
80 - 23 - 30
- 67 cm



Backrest angle
0 - 70°



Thigh angle
0 - 33°

Colours



Steel



Beech



Complies to the safety standard,
IEC 60601-2-52.



ISO 13485 certified company

Tel: +41 61 487 70 70
globalcomms@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
CH-4147 Aesch BL
SWITZERLAND