

Adeo



Our next generation of toileting solutions

Adjustable, supportive



Yes, you can.®

Introducing Adeo C407A a toilet frame

The Invacare Adeo C407A is a compact and lightweight toilet frame. With integrated adjustability the toilet frame can be adjusted to suit individual user needs. The easy to operate height-adjustable frame provides improved support whilst sat on the toilet.

Designed to be space-saving and fit into smaller toilet and bathroom environments whilst also being corrosion-free to ensure long term durability.



Key features:

- > Robust and supportive toilet frame
- > Easily adjustable for perfect support height
- > Compact design and fit into smaller toilet and bathroom environment
- > Provides a convenient and safer toileting alternative for users with reduced mobility.

Technical data

For more information about this product, including the user manual, please visit our website: www.invacare.eu



Adjustment steps
for height ▼

640 -740 mm



Total width ▼

590 mm



Total depth ▼

440 mm



Armrest length ▼

330 mm




Total product
weight ▼

2 kg



Maximum user
weight ▼

100 kg



Silver

Frame colour



Light grey

Plastic colour

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland



Yes, you can.®