

AQUATEC®

Dot



Needs little space and provides
great support



Yes, you can.®

The *Aquatec Dot* Shower Stool is the right choice if you have limited space in your bathroom and / or your shower. The shower stool convinces with its comfortable seat that provides cleaning access on both sides and also with its new and fresh design. The stool is very light weight and has a very easy to operate seating height adjustment. The rubber feet features a stable stand in- and outside the shower.

- Very easy to operate seating height adjustment in 7 steps from 40 cm to 55 cm
- Low product weight of 1,65 kg
- Tool-less assembly and very easy to clean
- Absolutely corrosion free - completely made out of plastic
- High load capacity: 135 kg



Features and Options



Modern shape

Comfortable seat with cleaning access on both sides.



Rubber feet

Ensure a stable stand in- and outside the shower.



Height adjustment

Intuitive operation.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

| | <i>Aquatec Dot</i> |
|--------------------------------------|--------------------|
| Shower Stool Single Pack | 1564679 |
| Shower Stool Pack of 2 - Retail Pack | 1564759 |
| Seat width | 343 mm |
| Seat depth | 213 mm |
| Seat height | 400 - 550 mm |
| Total width | 365 mm |
| Total depth | 365 mm |
| Total height | 400 - 550 mm |
| Load capacity | 135 kg |
| Weight | 1,65 kg |
| Colour | Grey |

Tel: +44 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



Yes, you can.®



© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU - 07/23