

H340 Izzo

advanced dual-function toileting solution



Next generation of toileting solutions

Dual-function, adjustable, supportive



Yes, you can.®

Introducing H340 Izzo an advanced dual-function toileting solution

The Invacare H340 Izzo can be used over the toilet as toilet frame or toilet seat raiser as well as a commode with the toilet pan with lid – depending on your individual requirements.

The easy to operate flip down armrests and the height adjustable seat support the transfer.










Key features:

- > Easy to handle seating height adjustment in 5 steps from 45 cm to 55 cm
- > Complete with elongated toilet seat with lid and pan with lid
- > Tool-less assembly
- > High load capacity: 113 kg

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

					
	Seat width ▼	Seat depth ▼	Seat height ▼	Total width ▼	Frame colour
H340 Izzo	36 cm	43 cm	45-55 cm	47 cm	White
					
	Total depth ▼	Toilet height ▼	Total product weight ▼	Maximum user weight ▼	Plastic colour
H340 Izzo	53 cm	82 - 92 cm	8,8 kg	113 kg	Grey

Tel: +44 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



Yes, you can.®