

H440 Omega Eco



Yes, you can.®

The *Invacare H440 Omega Eco* is the right choice if you want to combine function, comfort and design. With the snap button height adjustment and the innovative toilet recess insert, this toilet chair is simple and effortless to use for everyone.

- Modern and light design for the perfect fit in the living environment
- Adjustable seating height from 405 mm to 585 mm – simple to operate
- Easy handling with the recess insert
- Perfect over the toilet position
- Easy to clean cover material
- High load capacity: 140 kg

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Order numbers	H440 Omega Eco
Commode chair, 5 pack	1542903
Dimensions	
Seat width	450 mm
Seat depth	425 mm
Seat height	405 - 585 mm
Total width	585 mm
Total depth	515 mm
Total height	790 - 970 mm
Height to armrests	620 - 800 mm
Height, ground to lower seat edge:	
With holder for toilet pan	335 - 515 mm
Without holder for toilet pan	345 - 525 mm
Height floor to seat plate	405 - 585 mm
Width between the armrests	470 mm
Width, cut out in seat plate	200 mm
Depth, cut out in seat plate	250 mm
More info	
Load capacity	140 kg
Weight (kg)	7,5 kg
Colour	Lavender blue

Tel: +44 61 487 70 70
 hqeuropa@invacare.com
 www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



Yes, you can.®