

# Platinum<sup>®</sup> 9



Reduce cost with Invacare high  
flow concentrator



Yes, you can.™

The *Invacare Platinum 9* offers providers a value added and reliable global solution to treat patients with high flow requirements. Use as a stand alone device or couple together with *HomeFill* for complete patient care. While optimising oxygen concentration, *Platinum 9* minimises costs associated to maintenance and delivery. It allows patients to stay at home comfortably and safely thanks to a series of integrated alarms. Users will appreciate its simplicity of use combined with maximum security.



## Features and options

### High performance and reliability combined with simplicity

Powerful and reliable, *Invacare Platinum 9* is also user friendly. With its 9 litre oxygen capacity and 96-87% purity, *Platinum 9* is dedicated to patients with high oxygen needs. Easy to manoeuvre thanks to its integrated handle and casters, *Platinum 9* is also easy to use and guarantees extra safety for its user.



#### Filter

Filter is accessible without tools from the concentrator side thanks to door easy to open.



#### Homefill connection port

*Platinum 9* is compatible with *Homefill* system for higher flow above 5 LPM.





### Value added associated to cost control

Cost linked to inventory is kept to a minimum thanks to the small number of spare parts and specific accessories. With no maintenance for three years, Invacare offers you the best value added concentrator.

## Technical data




For more comprehensive pre-sales information about this product, including the product's user manual, please refer to your local Invacare website.

	Length	Width	Height	Weight	Charger	Power consumption	Pulse dose delivery	Outlet pressure
<i>Platinum 9</i>	 365 mm	 467 mm	 670 mm	 24 kg	 230 V, 50 Hz	 475 W	 1 - 9 l/min	 62 kPa

	Sound level	O2 purity	Altitude	ISO Compliant
<i>Platinum 9</i>	 < 47 dB(A)	 87-96%	 Up to 1230m	 ISO80601-2-69:2014

### SensO<sub>2</sub>-System

### Service interval

	O <sub>2</sub> -concentration
	higher than 85% - green light - <b>Normal</b>
	between 73% and 85% - yellow light - <b>Check filter and flow rate</b>
	lower than 73% - red light & audible alarm - <b>Call service technician</b>

No provider routine preventative maintenance for 3 years, replace filters as required.



Scan here to access the website  
<http://oxygen-therapy.net/en/>

Tel: +41 61 487 70 70  
hqeurope@invacare.com  
www.invacare.eu.com

Invacare International GmbH  
Benkenstrasse 260  
4108 Witterswil  
Switzerland



Yes, you can.