

AQUATEC®

# Pure Bidet



Modern, luxurious and hygienic



Yes, you can.®

The *Aquatec Pure Bidet* is ideal for those who experience difficulty with **personal hygiene**, or for those simply seeking **added luxury** at home. The modern looking bidet **fits easily** and **discreetly** to existing toilets and provides a gentle yet effective **cleansing solution** using a warm water spray. Key features include:

- **High functionality**
- **Inbuilt technology**
- **Fits to existing toilet**

### High functionality

- **Intuitive** hand control with large buttons
- **Adjustable** water temperature with **aerated** spray and variable jet settings
- **Personal** spray head positioning for both posterior and feminine settings
- Adjustable **warm air dryer**
- Temperature controlled **heated seat**
- **LED** night light
- **Soft close** seat and lid



### Inbuilt technology

- **Continuous** water heating for **instant** warm water that will not run out during cleaning
- Automatic spray head cleansing before and after use to **maintain cleanliness** of the system
- Equipped with a **thermal protector** to stop overheating
- **Automatic shut down** for energy saving which is activated via its seat sensor
- **Self-diagnosis** feature to automatically detect malfunction
- Earth Leakage Breaker (ELB) for **added safety**, giving peace of mind



### Fits to existing toilet

- Uses **standard pan fixing points**
- **Quick** and **easy** to assemble and detach
- WRAS and DIN EN 1717 conformity to ensure **safe connection** to public water pipes



### Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Total width	389 mm	Type of protection, entire device	IP X4
Total length	515 mm	Battery type remote control	2 x AA (LR6) / 1.5 V
Total depth	156 mm	Water supply pressure	79 - 735 kPa
Product weight	3.9 kg	Water supply temperature	3 - 35 °C
Max. user weight	125 kg	Water flow rate CLEANSING	0.7 l/min
Length of power cord	1.4 m	Water flow rate BIDET	0.85 l/min
Power supply	220 - 240 V AC, 50/60 Hz	Water pressure adjustment	3 levels
Max. power consumption	1725 - 2053 W		

Tel: +44 61 487 70 70  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
 Neuhofweg 51  
 4147 Aesch BL  
 Switzerland



© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU - 07/23

**Yes, you can.®**