


# Portable Ramps

Deployed in seconds



Yes, you can.®

# Portable Ramps

Invacare Portable Ramps are manufactured with a unique combination of uni-directional glass and carbon fibre materials to create an extremely strong and lightweight ramp suitable for a wide range of uses.

- > **Portable** - Quick and easy to fold into a single compact unit for easier storage or transportation in any vehicle or home environment.
- > **Lightweight** - Made from a unique combination of uni-directional glass and carbon fibre materials to give a strong, yet lightweight ramp that is easy to handle.
- > **Compact** - Multi-paneled design with a patented hinge system ensure the ramps fold in a concertina action giving a smaller width when folded.
- > **Safety** - To help improve safety during wheelchair transfer, our ramps are designed with pads at the top and bottom of the base offering a non-slip surface and feature clearly marked grip tape for improved visibility.
- > **Strong** - Tested to AS/NZ Standards with up to 300 kg weight load for peace of mind.
- > **Simple to use** - No assembly required, allowing quick and easy access for immediate use.



Invacare Portable Ramps are available in 4 models, in a wide range of sizes;

- > **Personal** - Portable, lightweight ramp that allows easier access in and around the home or workplace. Ideal for climbing kerbs or similar obstacles.
- > **Multipurpose** - Longer and wider, designed to cope with a wider range of obstacles making it more suitable for commercial and travel situations. Invaluable when venturing outdoors.
- > **Senior** - Ideally used where the doorway is narrow and a great solution for loading unoccupied equipment.
- > **Edge Barrier Limiter (EBL)** - Ideal for people who have busy lifestyles and require frequent support accessing, shops, places of work and other indoor environments. A 5cm Edge Barrier helps improve safety and security, helping to prevent any veering off course which is especially important with longer length ramps.
- > **Edge Barrier Limiter Xtra** - Designed with a 5cm Edge Barrier for improved safety and security, helping to prevent steering off course. The dual panel design folds in a bi-fold action. Up to 400kg weight load.



- > **Carry bag** for Personal and one of the Multipurpose models.



- > Marked 'TOP' to show head of ramp, grip tape, pads at both ends.



- > **Multi-paneled, hinge system** folds with a concertina action, allowing it to be folded to a single compact unit.

## Technical data ▼

For more information about this product, including the user manual, please visit our website: [www.invacare.com](http://www.invacare.com)

| Model/ Length cm  | Width  | Weight | Model/ Length cm | Width | Weight | Model/ Length cm | Width | Weight |
|-------------------|--------|--------|------------------|-------|--------|------------------|-------|--------|
| Personal 070*     | 74.5cm | 4.7kg  | Senior 090       | 75cm  | 4.8kg  | EBL 090          | 82cm  | 6kg    |
| Multipurpose 087* | 83.5cm | 5.6kg  | Senior 120       | 75cm  | 8.0kg  | EBL 120          | 82cm  | 9.7kg  |
| Multipurpose 100  | 83.5cm | 6.5 kg | Senior 165       | 75cm  | 10.8kg | EBL 165          | 82cm  | 13kg   |
|                   |        |        | Senior 200       | 75cm  | 11.7kg | EBL 200          | 82cm  | 13.7kg |
|                   |        |        | Senior 300       | 75cm  | 21.1kg | EBL 300          | 82cm  | 27.3kg |
|                   |        |        | Senior 350       | 75cm  | 24.7kg | EBL 350          | 82cm  | 31.5kg |

\*: sold with a carry bag

Tel: +41 61 487 70 70  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
 Neuhofweg 51  
 4147 Aesch BL  
 Switzerland



Yes, you can.®