

AQUATEC®

Sorrento



Shower and Bathtub Collection



Yes, you can.™

The outstanding modular and modern *Aquatec Sorrento* range is a complete Shower and Bathtub collection, the Shower Chair, Stools as well as the Bathtub Swivel Seat and the Bathtub Seat include great function with great style and design. The non-slipping soft touch surfaces and the soft backrest provide a secure and very comfortable seat for the user.

The wide adjustment range enables an easy fit to every user. The all-corrosion free, modular shower chair range is the best choice for those who need assistance in the bathroom for their personal hygiene.

Non-slipping Soft Touch Surface

The *Aquatec Sorrento* is the first shower chair with the innovative TPE Soft Touch Surface on the seat and the armrests. Not only the great look but also the comfortable feeling will persuade the user.

Adjustable soft backrest

The soft and adjustable backrest of the *Aquatec Sorrento* provides a very comfortable feeling for the user. With the wide adjustment range every user will feel pleased due to an easy personal adaption.



Sorrento as Shower Chair



Sorrento as Shower Stool



Sorrento as Bathtub Swivel Seat

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	<i>Sorrento Shower Chair</i>	<i>Sorrento Shower Stool w/armrests</i>	<i>Sorrento Shower Stool</i>	<i>Sorrento Shower Chair w/o recess</i>	<i>Sorrento Bathtub Swivel Seat</i>
Shower Chair	1532633			1532730	
Shower Stool w/armrests		1532660			
Shower Stool			1532661		
Bathtub Swivel Seat					1532693
Measurements					
Seat width (mm)	400	400	400	400	400
Seat depth (mm)	400	400	400	400	400
Seat height (mm)	460 - 610	460 - 610	460 - 610	460 - 610	
Total width (mm)	520	520	520	520	710 - 760
Total depth (mm)	430	430	430	430	650
Total height (mm)	850 - 990	655-805	460-610	850 - 990	520
Width between armrests (mm)	450	450		450	450
Armrest length (mm)	280	280		280	280
Armrest height (mm)	195	195		195	195
Backrest width (mm)	500			500	500
Backrest height (mm)	380			380	380
Backrest cover height (mm)	280			280	280
Width hygiene recess (mm)	140	140	140		140
Depth hygiene recess (mm)	170	170	170		170
High load capacity	135 kg	135 kg	135 kg	135 kg	130 kg
Weight	5,5 kg	4 kg	3,5 kg	5,5 kg	6 kg
Colour					



Tel: +44 61 487 70 70
 hqeuropa@invacare.com
 www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU - 07/23

Yes, you can.®