

Action^{®2}
NG



The Next Generation has arrived!



Yes, you can.

The **Action²_{NG}** was developed from the tried-and-tested **Action³** concept. The specifications meet the basic needs of the users whilst refusing to compromise on safety or comfort. The **Action²_{NG}** shares the same essential platform as the **Action³**, offering a low total weight for easy handling combined with excellent driving abilities. Due to the shared use of common components, the **Action²_{NG}** perfectly fits into the **Action_{NG}** Family.

Light and comfortable

With a weight from 16,2 kg, the *Invacare Action²_{NG}* fulfills all essential functional requirements. The stable folding frame means the wheelchair is easy to handle and is light weight. Cross-compatibility with other members of the **Action_{NG}** Family ensures that the **Action²_{NG}** can be easily adapted to meet individual user needs.

Quality of life

The tried-and-tested technology of the **Action_{NG}** Family products allow the users to go wherever they desire. Cross-compatibility with **Action³_{NG}** and **Action⁴_{NG}** make this an ideal, basic, fleet product.

Features and Options



Castor housing

The castor housing is mounted outside to have more space for the feet when transferring of foot propelling.



Adjustable center of gravity

The newly-designed rear wheel plate can be adjusted by 50 mm in depth and propose 3 positions in height.



Anti-tippers

Provide extra safety on slopes or uneven ground.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

| | Seat width | Seat depth | Seat angle | Seat height | Backrest angle | Backrest height | Front frame angle | Length of footrests |
|------------------------|--------------------------------|--------------|------------|--------------------|----------------|-----------------|-------------------|---------------------|
| Self Propulsion | 380 / 405 / 430 / 455 / 480 mm | 400 / 430 mm | 3° | 460 / 485 / 510 mm | 7° / 0° - 30° | 430 mm | 80° / 0° - 80° | 330 / 460 mm |
| Transit | 380 / 405 / 430 / 455 / 480 mm | 400 / 430 mm | 3° | 485 / 510 mm | 7° / 0° - 30° | 430 mm | 80° / 0° - 80° | 330 / 460 mm |

| | Armrest height | Driving unit width | Total width, folded | Total length inc. Footrests | Total length without footrests | Total height | Max. Safe slope | Radius |
|------------------------|------------------------------------|---------------------|---------------------|-----------------------------|--------------------------------|--------------------|-----------------|--------|
| Self Propulsion | 200 - 280 mm (n°1) 240 mm (n°2) | Seat width + 200 mm | 325 mm | 1030 / 1080 mm | 820 / 870 mm | 920 / 945 / 970 mm | 9° | 870 mm |
| Transit | 200 - 280 mm (n°1) 240 mm (n°2) | Seat width + 145 mm | 240 mm | 950 mm | 725 mm | 945 / 970 mm | 9° | 870 mm |

| | Max. User weight | Total weight | Mass of the heaviest part |
|------------------------|------------------|--------------|---------------------------|
| Self Propulsion | 125 kg | 16,2 kg | 8,9 kg |
| Transit | 125 kg | 14,9 kg | 12,1 kg |

Frame colour



Upholstery

Please notice colors may vary slightly from those displayed above.

Tel: +44 61 487 70 70
hqueurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland

