

Action^{®2}
NG



The Next Generation has arrived!



Yes, you can.

The **Action²_{NG}** was developed from the tried-and-tested **Action³** concept. The specifications meet the basic needs of the users whilst refusing to compromise on safety or comfort. The **Action²_{NG}** shares the same essential platform as the **Action³**, offering a low total weight for easy handling combined with excellent driving abilities. Due to the shared use of common components, the **Action²_{NG}** perfectly fits into the **Action_{NG}** Family.

Light and comfortable

With a weight from 16,2 kg, the *Invacare Action²_{NG}* fulfills all essential functional requirements. The stable folding frame means the wheelchair is easy to handle and is light weight. Cross-compatibility with other members of the **Action_{NG}** Family ensures that the **Action²_{NG}** can be easily adapted to meet individual user needs.

Quality of life

The tried-and-tested technology of the **Action_{NG}** Family products allow the users to go wherever they desire. Cross-compatibility with **Action³_{NG}** and **Action⁴_{NG}** make this an ideal, basic, fleet product.

Features and Options



Castor housing

The castor housing is mounted outside to have more space for the feet when transferring of foot propelling.



Adjustable center of gravity

The newly-designed rear wheel plate can be adjusted by 50 mm in depth and propose 3 positions in height.



Anti-tippers

Provide extra safety on slopes or uneven ground.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat width	Seat depth	Seat angle	Seat height	Backrest angle	Backrest height	Front frame angle	Length of footrests
Self Propulsion	380 / 405 / 430 / 455 / 480 mm	400 / 430 mm	3°	460 / 485 / 510 mm	7° / 0° - 30°	430 mm	80° / 0° - 80°	330 / 460 mm
Transit	380 / 405 / 430 / 455 / 480 mm	400 / 430 mm	3°	485 / 510 mm	7° / 0° - 30°	430 mm	80° / 0° - 80°	330 / 460 mm

	Armrest height	Driving unit width	Total width, folded	Total length inc. Footrests	Total length without footrests	Total height	Max. Safe slope	Radius
Self Propulsion	200 - 280 mm (n°1) 240 mm (n°2)	Seat width + 200 mm	325 mm	1030 / 1080 mm	820 / 870 mm	920 / 945 / 970 mm	9°	870 mm
Transit	200 - 280 mm (n°1) 240 mm (n°2)	Seat width + 145 mm	240 mm	950 mm	725 mm	945 / 970 mm	9°	870 mm

	Max. User weight	Total weight	Mass of the heaviest part
Self Propulsion	125 kg	16,2 kg	8,9 kg
Transit	125 kg	14,9 kg	12,1 kg

Frame colour



Upholstery

Please notice colors may vary slightly from those displayed above.

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland

