

# Action<sup>®1</sup> R



A feature rich, entry level chair



Yes, you can.™

The **Action<sup>1</sup> R** is an entry level wheelchair that doesn't compromise on functionality. Built using decades of in-house design experience, the **Action<sup>1</sup> R** will be available in three versions; a self-propelling with 24" spoked rearwheels and aluminium handrims version, a self-propelling 24" rearwheels version with drumbrakes and a transit version with small 12" rearwheels. With

six standard sizes to choose from, it is quick and easy to find the right size. The removable, swing away armrests are easy to remove and the armpad can be used in two positions. The **Action<sup>1</sup> R** offers:.

- **Improved function and comfort**
- **Six sizes and three versions**
- **Proven quality from a market leader**

## Improved functionality

New armrests swing out the way for transfers and offer two armpad positions - desk and full length - which can be removed quickly and easily for transport.



### Proven reliability and quality

With years of experience and using tried-and-tested technology, the **Action<sup>1</sup> R** is a reliable partner, the one you can count on.



### Elevating legrests

The optional Eco elevating legrests provide **extra comfort for users**.



### Cane holder

This optional cane holder can be mounted to securely hold a cane or crutches.



### Anti-tippers

Dynamic anti-tippers provide **extra safety on slopes** or when backrest is reclined.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat width	Seat depth	Seat angle	Seat height	Backrest angle	Backrest height	Front frame angle	Length of footrests
Self Propulsion	380 / 405 / 430 / 455 / 480 / 505 mm	425 mm	3°	495 mm	0° / 7°	430 mm	70°	330 - 455 mm
Transit	380 / 405 / 430 / 455 / 480 / 505 mm	425 mm	3°	495 mm	0° / 7°	430 mm	70°	330 - 455 mm

	Armrest height	Driving unit width <sup>4</sup>	Total width, folded	Total length inc. footrests <sup>4</sup>	Total length without footrests	Total height	Max. safe slope <sup>1</sup>	Radius <sup>2</sup>
Self Propulsion	240 mm	608 - 700 mm (Seat width +195 mm)	250 mm	1085 - 1125 mm	835 - 875 mm	950 mm		
Transit	240 mm	578 - 670 mm (seat width +160 mm)	260 mm	955 mm	710 mm	950 mm		

	Max. user weight	Total weight	Mass of the heaviest part <sup>3</sup>
Self Propulsion	125 kg	from 17.5 kg	from 10 kg
Transit	125 kg	from 17.5 kg	from 10 kg

**1. Stability Values**  
 Max. slope parking brake: 7°  
 Static stability downhill: 18°  
 Static stability uphill: 9°  
 Static stability sideways: 19°

**2. Radius**  
 Mini turning radius: 915 mm  
 Pivot width: 1280 mm


**3. Folded Dimensions**  
 Folded height: 735 - 950 mm  
 Folded depth: 1025 mm  
 Folded width: 250 mm

**4. Overall dimensions**  
 In some configurations, overall width or length exceed the authorised limits, access to emergency escape routes is not possible




This wheelchair is not intended to be used as a seat in a motor vehicle

## Frame colour



## Upholstery



Please notice colors may vary slightly from those displayed above.

Tel: +44 61 487 70 70  
 hqeuropa@invacare.com  
 www.invacare.eu.com

Invacare International GmbH  
 Neuhofweg 51  
 4147 Aesch BL  
 Switzerland



Yes, you can.®