

# Alu Lite



A lighter way to get around



Yes, you can.

The *Invacare Alu Lite* is the basic-entry wheelchair in the Invacare manual range. The lightweight aluminium construction of this transit wheelchair makes it easy to manoeuvre, even for an elderly carer, as well as providing a wheelchair that is easily transportable.

### Easy to transport and store

Being lightweight (aluminium chassis), compact and with a foldable backrest, the *Invacare Alu Lite* is very easy to handle, transport and store. Simply remove the footrests, fold the backrest and chassis and store away.

### Pushing? No problem!

The overall manoeuvrability of a wheelchair is crucial, which is why the *Alu Lite* is designed to be very easy to push and manoeuvre, both indoor and outdoor. Due to the smaller overall dimensions and general compactness of this model, getting around in tight spaces is never a problem.

### Features and Options



#### Folding backrest

Reduces overall size when folded for transportation or storage.



#### Swing away footrest

Footplates fold up as shown and footrests swing away, allowing easier access in or out of the chair.



#### Foldable cross frame

Lightweight frame folds flat for easier portability.

### Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.


	Seat width	Seat depth	Seat angle	Seat height	Backrest angle	Backrest height	Front frame angle	Length of footrests
<i>Invacare Alu Lite</i>	405 / 455 mm	430 mm	0°	490 mm	0°	450 mm	80°	340 - 460 mm

	Armrest height	Driving unit width	Total width, folded	Total length inc. Footrests	Total length without footrests	Total height	Max. Safe slope	Radius
<i>Invacare Alu Lite</i>	230 mm	Seat width + 110 mm	250 mm	1030 mm	760 mm	945 mm	12,5°	890 mm

	Max. User weight	Total weight	Mass of the heaviest part
<i>Invacare Alu Lite</i>	100 kg	12 kg	10,2 kg


### Frame colour

Blue



### Upholstery

Black nylon



Please notice colors may vary slightly from those displayed above.

Tel: +44 61 487 70 70  
 hqeuropa@invacare.com  
 www.invacare.eu.com

Invacare International GmbH  
 Neuhofweg 51  
 4147 Aesch BL  
 Switzerland



**Yes, you can.®**