



Invacare[®] | **Mattress**

Therapeutically designed to
give a great night's sleep



The design of this mattress allows even weight distribution by conforming to the shape and form of the body, promoting comfort and support.

Air is allowed to circulate within the mattress, minimising the build-up of heat and moisture and delivering a comfortable night's sleep.

Available in two cover options.

White fabric cover option


Luxurious white fabric cover for softness and durability.

Blue fabric cover option

Blue two-way stretch breathable water resistant polyurethane fabric.



Features and Options



Double sided – regular turning will extend the life of the mattress.



Castellated surface for comfort and therapeutic support.



Optional water resistant liner available.








White fabric cover.



Blue water resistant fabric.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Mattress sizes	Height	Total product weight	Max. user weight	Cleaning
				
2000 x 855 mm and 1950 x 880 mm	140 mm	13 kg	140 kg	Blue cover: 80° White fabric Cover: 40°

2 year warranty


Mattress is manufactured to:

- European Standard: BS EN 597-1 Cig
- European Standard: BS EN 597-2 Match
- BS7177:1996 Crib 5. (Blue cover only)

Frame colours

White fabric

Baltic blue



QR Code

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland

