

# Bedco Backrest



Simple solution – big impact



Yes, you can.

This electrically operated backrest from Invacare is designed for use in private homes in standard private beds. Easy to use and activate, the overlay mattress stays in place thanks to the fabric's texture.

The **Bedco Backrest** is ideal for travelling with, as it is very easy to mount and remove.

## Features and Options



### Motor can be mounted on either side

No loose parts: the motor, handset and transformer are always connected. The motor can easily be mounted on both sides.



### Bag for easy transport

A handy bag for easy storage of hand set and transformer.



### Bedco Anti-Sliding

Prevents the backrest from sliding when the backrest is being operated.

### Comfort and functionality

The *Invacare Bedco Backrest* comes with maximum functionality. The backrest helps the client to establish a sitting position – making it easier to get in and out of bed.

The texture of the fabric ensures that the mattress stays firmly in place. The backrest is available in 4 widths to meet individual needs.

### Suitable for transport

Easy to transport (light weight, and no loose parts) the backrest can easily be mounted without need of additional tools.

The motor can be placed on either side to meet individual needs.

A 24V system increases the levels of safety.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Width outside / inside	Backrest angle	Total weight	Max. user weight	Voltage	Max. cleaning temperature
750 mm	800 / 700 mm	70°	9,5 kg	135 kg	24V	60°
900 mm	950 / 850 mm	70°	10 kg	135 kg	24V	60°
1050 mm	1100 / 1000 mm	70°	10,5 kg	135 kg	24V	60°
1200 mm	1450 / 1200 mm	70°	13 kg	150 kg	24V	60°

## Frame colour

Aluminium grey

Khaki



Tel: +41 61 487 70 70  
 hqeuropa@invacare.com  
 www.invacare.eu.com

Invacare International GmbH  
 Neuhofweg 51  
 4147 Aesch BL  
 Switzerland



Yes, you can.®



© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU - 04/23