

AQUATEC®

Dot



Needs little space and provides  
great support



Yes, you can.®

The *Aquatec Dot* Shower Stool is the right choice if you have limited space in your bathroom and / or your shower. The shower stool convinces with its comfortable seat that provides cleaning access on both sides and also with its new and fresh design. The stool is very light weight and has a very easy to operate seating height adjustment. The rubber feet features a stable stand in- and outside the shower.

- Very easy to operate seating height adjustment in 7 steps from 40 cm to 55 cm
- Low product weight of 1,65 kg
- Tool-less assembly and very easy to clean
- Absolutely corrosion free - completely made out of plastic
- High load capacity: 135 kg



## Features and Options



### Modern shape

Comfortable seat with cleaning access on both sides.



### Rubber feet

Ensure a stable stand in- and outside the shower.



### Height adjustment

Intuitive operation.

## Technical data

**For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.**

	<i>Aquatec Dot</i>
Shower Stool Single Pack	1564679
Shower Stool Pack of 2 - Retail Pack	1564759
Seat width	343 mm
Seat depth	213 mm
Seat height	400 - 550 mm
Total width	365 mm
Total depth	365 mm
Total height	400 - 550 mm
Load capacity	135 kg
Weight	1,65 kg
Colour	Grey

**Tel: +41 61 487 70 70**  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
 Neuhofweg 51  
 4147 Aesch BL  
 Switzerland



**Yes, you can.®**



© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU - 07/23