


Essential Basic

Quality foam mattress

A comfortable and supportive surface
for people who do not require
pressure injury prevention.



A cost-effective solution

For comfort and support



Yes, you can.®

Affordable quality mattress



Comfort & Support

- › Durable quality foam provides comfortable support.
- › One-piece foam allows the mattress to be turned inside the cover for increased product longevity.
- › Ideal for home, residential and long-term care environments.



Multi-Directional Stretch

- › Multi-stretch, water-resistant, polyurethane fabric cover helps reduce shear and friction forces.
- › Zipped around 3 sides for ease of inspection and laundering.
- › White inner liner allowing easy inspection for strike-through or other cover damage.



Quality Assured

- › Quality foam for durability.
- › Excellent recovery properties.
- › Manufactured to fire standards BS EN 597-1 Cig, BS EN 597-2 Match and BS 7177: Crib Source 5.

Technical data ▼

For more information about this product, including the user manual, please visit our website: www.invacare.com



Sizes available ▼

2000 × 880 × 150 mm
2050 × 880 × 150 mm



Max. User Weight ▼

140 kg



Max. Cleaning Temperature ▼

95°

Manufactured to: BS EN 597-1 Cig, BS EN 597-2 Match and BS EN 7177: Crib Source 5

Essential nursing care is pivotal in pressure injury management. This mattress may positively contribute to the outcome of a pressure injury prevention plan that includes education, clinical judgment and action planning based on vulnerability, which are fundamental factors in the prevention of pressure injuries. A range of assessment scales can be used as a formal method of assessing risk from pressure injury development, and should be used in conjunction with an informal assessment (informed nursing judgement) as part of a comprehensive and individualised plan.

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland



Yes, you can.®



© 2024 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Essential Basic - EU - 01/2024