

Essential Care

Quality castellated foam mattress

A comfortable pressure-redistributing support surface for people at risk of developing pressure injuries.



A cost-effective solution

For comfort and pressure redistribution



Yes, you can.®

Affordable quality mattress



Pressure Redistribution

- › Castellations help redistribute weight over a greater area.
- › Durable castellated foam for immersion and envelopment to provide optimal pressure redistribution.
- › Excellent recovery properties ensuring good pressure redistribution.
- › Suitable for those considered 'At Risk' to 'High Risk' of developing pressure injuries.



Shear & Friction Reduction

- › Horizontal castellations help reduce shear forces when moving up and down the bed.
- › Multi-stretch, water resistant, polyurethane fabric cover helps further reduce shear and friction forces.
- › Smooth cover supports in the reduction of frictional forces.



Comfort

- › Castellations designed to improve comfort and support.
- › Three-sided zipped cover with white inner liner for easy strikethrough inspection and removal for laundering.
- › Suitable for multiple care environments.

Technical data ▼

For more information about this product, including the user manual, please visit our website: www.invacare.co.uk



Sizes available ▼



Max. User Weight ▼



Max. Cleaning Temperature ▼

150 mm height

1950 × 880 × 150 mm
2000 × 880 × 150 mm
2050 × 880 × 150 mm

140 kg

95°

Manufactured to: BS EN 597-1 Cig, BS EN 597-2 Match and BS EN 7177: Crib Source 5

Essential nursing care is pivotal in pressure injury management. This mattress may positively contribute to the outcome of a pressure injury prevention plan that includes education, clinical judgment and action planning based on vulnerability, which are fundamental factors in the prevention of pressure injuries. A range of assessment scales can be used as a formal method of assessing risk from pressure injury development, and should be used in conjunction with an informal assessment (informed nursing judgement) as part of a comprehensive and individualised plan.



T: +44 1656 776222
uk@invacare.com
www.invacare.co.uk
www.invacare.ie

Invacare UK & Ireland
Unit 4 Pencoed Technology Park
Bridgend CF35 5AQ



Yes, you can.®



© 2024 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Essential Care - UK - 01/2024