


Essential Visco

Quality dual layer foam mattress

A pressure redistributing support surface providing additional immersion and envelopment.



A cost-effective solution

For pressure redistribution



Yes, you can.®

Affordable quality mattress



Pressure Redistribution

- › Viscoelastic foam top layer conforms to the body shape providing additional immersion and envelopment.
- › Quality foam base delivers increased support and reduces the likelihood of bottoming out.
- › Excellent recovery properties ensuring good pressure redistribution.



Shear & Friction Reduction

- › Multi-stretch, water-resistant, polyurethane fabric cover helps reduce shear and friction.
- › White inner liner allowing easy inspection for strikethrough or other cover damage.
- › Smooth cover supports in the reduction of frictional forces.



Quality Assured

- › Quality dual layer foam for durability and comfort.
- › Manufactured to fire standards BS EN 597-1 Cig, BS EN 597-2 Match and BS 7177: Crib Source 5.
- › Suitable for all care environments.

Technical data ▼

For more information about this product, including the user manual, please visit our website: www.invacare.com



Sizes available ▼



Max. User Weight ▼



Max. Cleaning Temperature ▼

	Sizes available ▼	Max. User Weight ▼	Max. Cleaning Temperature ▼
150 mm height	1950 x 880 x 150 mm	140 kg	95°
120 mm height	1950 x 880 x 120 mm	110 kg	95°

Manufactured to: BS EN 597-1 Cig, BS EN 597-2 Match and BS EN 7177: Crib Source 5

Essential nursing care is pivotal in pressure injury management. This mattress may positively contribute to the outcome of a pressure injury prevention plan that includes education, clinical judgment and action planning based on vulnerability, which are fundamental factors in the prevention of pressure injuries. A range of assessment scales can be used as a formal method of assessing risk from pressure injury development, and should be used in conjunction with an informal assessment (informed nursing judgement) as part of a comprehensive and individualised plan.

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



Yes, you can.®



© 2024 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Essential Visco - EU - 01/2024