

# **Essential** Care

*Quality Foam Mattress*




A cost effective solution for comfort  
and pressure redistribution



The new **Essential Care** mattress is a pressure redistributing support surface for the community, residential and community healthcare environment. Its castellated surface provides a cost effective solution for comfort and pressure redistribution to meet the needs of individuals at risk or high risk of developing pressure ulcers.

- ▶ Cut foam castellations
- ▶ Multi-stretch polyurethane fabric cover
- ▶ Zipped around 3 sides for ease of inspection and laundering



### Cover

- ▶ The outer cover is a multi stretch, water resistant, polyurethane fabric providing excellent recovery and durability, which helps to reduce shear and friction forces. The cover also has a white inner liner allowing easy inspection for strikethrough or other cover damage.

### Cut out castellations

- ▶ Castellations improve comfort and ensure effective pressure redistribution.

### Design

- ▶ The **Essential Care** mattress has an innovative and modern cover design with the simple wave giving the mattress a fresh and current feel.

## Technical data ▼

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.



Size range ▼



Max. user weight ▼



Max. cleaning temperature ▼

Nominal density range (kg/m<sup>3</sup>) ▼

Article no. ▼

	Size range ▼	Max. user weight ▼	Max. cleaning temperature ▼	Nominal density range (kg/m <sup>3</sup> ) ▼	Article no. ▼
<b>Essential Care</b> (150 mm height)	2000 x 880 x 150 mm	140 kg	95°	32 - 34	1588619
<b>Essential Care</b> (120 mm height)	2000 x 880 x 120 mm	110 kg	95°	32 - 34	1588614

Mattresses is manufactured to:  
EN 597-1 Cig, BS EN 597-2 Match and BS 7177: Crib Source 5

1. *Essential nursing care is pivotal in pressure ulcer prevention. This mattress will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement).*

Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS -  
County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)



Yes, you can.®

