



Flexibility and comfort – an ideal solution for Acute and Residential care

The *Invacare Jasmine* Mobile Lifter is a flexible solution that offers a high level of comfort for all clients. A versatile lifter, the **Jasmine** provides safe patient handling up to 32 stone (200kg) weight capacity. With an expansive, wide-opening base and high lifting mechanism, the **Jasmine** offers easy manoeuvrability and positioning for a wide range of lifts. With the latest control box technology, the **Jasmine** offers a 'Soft Start' and 'Soft Stop' function that provides a more comfortable and dignified transfer.



200kg weight limit
5 year warranty

Trusted Solutions, Passionate People





Features and Options

Flexibility and optimized maneuverability

The *Invacare Jasmine* is exceptionally mobile and easy to transport. The wide-opening base allows easy access for a wheelchair or armchair, whilst the high lifting mechanism is essential for transfers to and from elevated surfaces. Flexible and with optimized maneuverability, the *Jasmine* covers the needs of both the carer and client.

Safety first, with an intelligent control box

The *Jasmine* is fitted with the control box with diagnostic system: an intelligent service lifetime monitoring system that highlights the battery status and service condition of the actuator to ensure safe patient handling.



Spreader bar options

The *Jasmine* can be ordered with 2 or 4-point spreader bars. In widths: 450 mm and 550 mm.



Scale

A scale is available. Approved according to medical standard EN 45501, certified and sealed.



Linak® Jumbo™ Care

The intelligent monitoring system includes a service reminder, emergency operation, audio battery alarm, and battery indicator.

Technical data



Invacare Jasmine 134 cm



69 cm



105 cm



54 - 192 cm



12 cm



Invacare Jasmine 153 cm



200 kg




43 kg



40 cycles

Frame colours



Moderate white



Blue dusk

Warranty

5 year warranty:
Frame, Control box,
Actuator

2 year warranty:
Hand control

1 year warranty:
Battery

Invacare Limited
Pencoed Technology Park -
Pencoed - Bridgend CF35 5AQ -
United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk
Invacare Ireland Ltd
Unit 5 - Seatown Business
Campus
Seatown Road - SWORDS -
County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie
© 2016 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Jasmine - UK - 01/2016
Brochure re-order code: 1539003



safepatienthandling



Yes, you can.®