

Kirilos and Atlante



*Invacare **Atlante** - electrical version*



*Invacare **Atlante** - manual version*

**Compact and easy to maneuver -
the best choice for homecare**



Yes, you can.®

Compact in size and offering maximum maneuverability, these mobile lifters provide the ideal solution for the homecare environment. Carers will appreciate the easy handling and smooth operating functions of the **Atlante** and **Kirilos** models.

The **Atlante** lift comes in both an electric and hydraulic version; **Kirilos** is purely electric. Both models come equipped with a handle for easy leg opening.


Features and Options

Easy handling



Spreader bar

Spreader bar rotates 360 degree for easy positioning of the patient.



For easy leg opening

Equipped with a handle for comfortable leg opening.

Universal Standard sling, size M, solid

The lifter comes with a Universal Std solid sling (size M). Additional slings are available with or without head support and with or without commode opening.

The **Atlante** and **Kirilos** lifters are designed to make life easier for carers when handling patients. The lightweight, compact nature of the lifters makes them easy to maneuver even when space is limited. The spreader bar rotates 360 degrees, providing maximum handling opportunities for the carer. In addition the lifters can be dismantled for easy transportation and storage.

High safety

The control box of the **Kirilos** and the electric version of **Atlante** are equipped with emergency lowering and emergency stop. The low battery status is displayed by an LED light as well as an acoustic signal.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

| | Overall length of base | Outside width: legs closed | Total product weight | Lifting range | Height of legs | Inside width: legs opened | Max. user weight | Number of cycles |
|-----------------------------------|------------------------|----------------------------|----------------------|---------------|----------------|---------------------------|------------------|------------------------|
| Invacare Atlante hydraulic | 1170 mm | 700 mm | 37 kg | 610 - 1530 mm | 150 mm | 870 mm | 140 kg | |
| Invacare Atlante | 1170 mm | 700 mm | 33 kg | 610 - 1530 mm | 150 mm | 870 mm | 140 kg | 40 cycles ¹ |
| Invacare Kirilos | 1200 mm | 680 mm | 33 kg | 575 - 1490 mm | 110 mm | 900 mm | 140 kg | 40 cycles ¹ |

| | Battery capacity | Turning diameter | Wheel diameter | Max. cleaning temperature |
|-----------------------------------|------------------|------------------|----------------|---------------------------|
| Invacare Atlante hydraulic | | 1230 mm | 100 mm | 95° C |
| Invacare Atlante | 2,9 Ah | 1230 mm | 100 mm | 95° C |
| Invacare Kirilos | 2,9 Ah | 1230 mm | 75 mm | 95° C |

Frame colour

Grey



1. without battery charge with batteries at 50% of full capacity

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland



Yes, you can.®