

AQUATEC®

# Kogia



Bathing simply reinvented



Yes, you can.®

The *Aquatec Kogia* makes bathing simple and has been designed with **comfort, ease of use** and **safety** in mind. Weighing only 10kgs, the **Kogia** is **lightweight** yet able to support up to 140kgs. Having been tested against the toughest safety standards, the bathlifter has a robust and stable structure allowing you to bathe with peace of mind.

- Easy to use with a smooth riser mechanism
- Easy to assemble and maintain
- Lightweight and easy to move and handle

## Features and Options



### Floating handset

This easy to use, floating hand control has been **ergonomically** designed and can be placed in its handset holder when not in use.



### Open base plate


There is no need to remove the **Kogia** from the bath when you clean thanks to the open base plate.

## Comfortable bathing experience

While straightforward in nature, the **Kogia** has many features including its easy to use hand-control, smooth moving riser mechanism, its ability to lower to within 60 mm of the base of the bath and its **shorter seat**, allowing more room to bathe.

## Assembly and maintenance

The **Kogia**'s simple design makes assembly and maintenance very easy. Cleaning is effortless due to its smooth surface and open frame base plate, allowing you to clean the **Kogia** while it is still in the bath. Its lightweight design makes it easy to move in and out of the bath, all of which allow for a **easy, comfortable** and **safe** bathing experience.



## Technical data

**For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.**

	<i>Aquatec Kogia</i>
Seat width with side flaps folded in	370 mm
Seat width with side flaps folded out	715 mm
Seat depth	450 mm
Seat height at highest position	420 mm
Seat height at lowest position	60 mm
Backrest width	370 mm
Height of backrest	660 mm
Base area	570 x 266 mm
Total weight without battery	10 kg
Max. load	140 kg
Backrest can be reclined by up to	40°
Battery	Lead storage battery

Tel: +41 61 487 70 70  
 hqeuropa@invacare.com  
 www.invacare.eu.com

Invacare International GmbH  
 Neuhofweg 51  
 4147 Aesch BL  
 Switzerland



**Yes, you can.®**

