

# H263/H273 Lima



The foldable shower commode that comes with everything



Yes, you can.®

The *Invacare Lima* is the perfect choice if you are looking for a space-saving and fully equipped shower commode. All accessories are already included and the easy-to-handle folding mechanism provides a perfect storage solution for institutional and home care demands.

The rust resistant aluminum framed commodes are ideal in the shower for users who need assistance for their personal hygiene.

### Space-saving foldable design

See how easy a shower commode can be folded. The *Invacare H263 Lima* and *Invacare H273 Lima* are ready for storage or transport in just a few seconds. This saves space in the bath- or repository room.

### Variable seating positions

The padded 4-position seat of the *Invacare Lima* provides maximum comfort for its user and can be used over every toilet. It can be easily adapted and perfect adjustment for individual patient needs is possible.



### Features and Options



#### Wheel options

The choice of 5" castors or 24" self propelling wheels provides a custom fit solution.



#### All corrosion free

The aluminum frame and stainless steel hardware makes it completely rust resistant.



#### Soft backrest

The fabric backrest shows maximum comfort to the patient.

### Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat height	Seat width	Total width	Total depth	Total height	Total width folded	Total weight	Max. user weight
<b>H263</b>	580 mm	460 mm	570 mm	1050 mm	1000 mm	225 mm	14 kg	135 kg
<b>H273</b>	580 mm	460 mm	680 mm	1120 mm	1000 mm	335 mm	18,5 kg	135 kg

### Colour

Blue



Tel: +41 61 487 70 70  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
 Neuhofweg 51  
 4147 Aesch BL  
 Switzerland



Yes, you can.®