

matrx[®]
SEATING SERIES

Flo-tech[®] Back Support



The Matrx Flo-tech Back Support
Range



The *Invacare* **Matrix Flo-tech Back Support** range has been designed to provide added stability to a wheelchair or standard armchair and it can be used in conjunction with any **Matrix Flo-tech** cushion.

The **Matrix Flo-tech Back Support** range provides excellent comfort and improved user positioning, which in turn helps to improve or maintain independence whilst in the seated position in a wheelchair or armchair. The unique supportive design offers protection from pressure ulcer formation and spinal deformities which can cause severe pain.

A highly protective Back Support The added support of the Lumbar

The **Back Support** provides a high degree of stability and protection from pressure ulcer development.

The **Back Support with Lumbar** can be used for those patients who wish to have increased lateral and pelvic stability. The Visco Elastic Lumbar Support contours to the individual and can be easily adjusted and secured with Hook&Loop pads, to accommodate a range of heights.

Features and Options



Cover is moisture resistant

Cover is two-way stretch, moisture-resistant and vapour-permeable, which helps to prevent a build up of shear and friction forces against the user.

Extendable Straps


Offers the user a securely fixed **Back Support** on a standard armchair or wheelchair. The straps are easy to adjust and can be fixed to a position that provides the user with the maximum comfort and support.

Ultra-fresh™ added to all Flo-tech foam for antimicrobial and odour protection

Antimicrobial treatment that has been added to all **Flo-tech** foam to help make it resistant to the growth of odour causing bacteria. Helps provide longer lasting freshness.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local *Invacare* website.



Matrix Flo-tech Back Support

165 kg


80°

1.1 kg

Back Support and Lumbar are manufactured to European Standards: BS EN 1021-1 Cig, BS EN 1021-2 Match and BS 5852: Crib 5

Colours

Black Smooth



Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland



Yes, you can.®