

Octave<sup>®</sup>



Setting new standards....



Yes, you can.<sup>®</sup>

With a safe working load up to 385 kg, the *Invacare Octave* is a state-of-the-art profiling carebed solution that assists in the care and comfort of bariatric clients. With a fully electrical 4-sectioned mattress support as standard, the **Octave** provides optimal comfort and support to the user. The **Octave** provides an ideal working environment for the carer. The construction of the mattress platform allows easier access to the client. The bed ends can be detached independently of the side rails for quick and easy access to both head and feet.

## Great comfort and design

*Invacare Octave* is designed specifically with the bariatric client in mind. The bed frame extends to the edges of the bed base for enhanced support.

An easy to use hand control allows the client to operate the bed height and profiling function. The **Octave** is a substantial bed, however its design includes all the well-known qualities of Scandinavian design.

## Features and Options



### Lifting pole

The handle can be adjusted in both height and depth.



### Extendable bedend

The bedend is extendable and offers flexibility. Choose between 10 cm to 20 cm extension.



### Support handle

Makes getting in and out of bed much easier for the patient.

## Delivery and storage

To make delivery and transportation as easy as possible, **Octave** is delivered on a transportation rack.

Dismounting and assembly of **Octave** couldn't be simpler despite the size of the bed. Furthermore the transportation rack makes it possible to economise on storage space and makes the **Octave** very easy to move between locations.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local **Invacare website**.

	Width outside	Length outside Length inside	Height adjustment	Mattress support dimensions	Backrest angle	Thigh angle	Legrest lift	Total weight
<i>Invacare Octave</i>	1190 mm	2210 mm 2070 mm Mattress: 2000 mm	400 - 800 mm Total height: 1720 - 2120mm	810-340- 330-520 mm (4-sectioned)	0 - 57°	0 - 30°	0 - 30°	Steel 180 kg Wood: 199 kg

	Max. user weight
<i>Invacare Octave</i>	SWL 385 kg 350 kg



Complies with IEC 60601-2-52


### NB:

The bed must not be used by patients under 12 years of age, or by patients with body size equivalent to an average 12 year old or smaller.

## Frame colour

Silver

Beech



Tel: +41 61 487 70 70  
hqeurope@invacare.com  
www.invacare.eu.com

Invacare International GmbH  
Neuhofweg 51  
4147 Aesch BL  
Switzerland



Yes, you can.®



© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU - 04/23