

# *Omega Commode*



Yes, you can.®

The *Invacare Omega* is the right choice if you want to combine function, comfort and design.



#### *Invacare Omega Eco H440*

- Stationary commode chair
- Space-saving, commode can be stacked
- Stable and corrosion-protected design with welded tubes
- Height adjustable seating height
- Perfect over the toilet position
- Easy handling with the recess insert

#### *Invacare Omega H450*

- Stationary commode chair
- Space-saving: commode chair can be stacked
- Stable and corrosion-protected design with welded tubes
- Easy handling: commode pan can be removed from the side and from above

#### *Invacare Omega H407*

- Foldable, Space saving commode
- Ideal for infrequent or travel use
- Contoured backrest
- Complete with sanitary pan

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Order numbers	H440 Omega Eco	Omega H450	Omega H407
Commode chair, 5 pack	1542903	1471163	-
Commode chair, 2 pack	-	-	1471010
Dimensions			
Seat width	450 mm	495 mm	440 mm
Seat depth	425 mm	455 mm	420 mm
Seat height	405 - 585 mm	500 mm	450 mm
Total width	585 mm	635 mm	580 mm
Total depth	515 mm	360 mm	560 mm
Total height	790 - 970 mm	890 mm	910 mm
Height to armrests	620 - 800 mm	690 mm	660 mm
Height, ground to lower seat edge:			
With holder for toilet pan	335 - 515 mm	370 mm	-
Without holder for toilet pan	345 - 525 mm	430 mm	-
Height floor to seat plate	405 - 585 mm	525 mm	475 mm
Width between the armrests	470 mm	510 mm	460 mm
Width, cut out in seat plate	200 mm	215 mm	215 mm
Depth, cut out in seat plate	250 mm	250 mm	255 mm
More info			
Load capacity	140 kg	130 kg	130 kg
Weight (kg)	7,5 kg	8,3 kg	8,5 kg
Colour	Lavender blue	Lavender blue	Lavender blue

Tel: +41 61 487 70 70  
 hqeuropa@invacare.com  
 www.invacare.eu.com

Invacare International GmbH  
 Neuhofweg 51  
 4147 Aesch BL  
 Switzerland

