

Invacare[®] | **Propad[™] Profile Cushion**

Profiled cut-foam cushion designed for comfort with castellations and support

Trusted Solutions,
Passionate People



Features and Options

The *Invacare Propad Profile cushion* provides individuals with an option of having a lower 3" high cushion, rather than a standard height (4") original cushion. Some users prefer a lower profile cushion, finding it more compatible with their daily living activities.



Cover helps to prevent a build up of shear and friction forces

Cover is two-way stretch, water-resistant and vapour permeable. Cover features an anti-slip base.

Available in a range of sizes

The *Invacare Propad Profile Cushion* is available in 102 sizes, ranging from 25 x 25 cm (250 x 250 mm) – 60 x 55 cm (600 x 550 mm).

Full length flap conceals the zip

The full length flap concealing the zip reduces the potential for fluid ingress into the cushion and protects the foam core. The cover can be removed and laundered.

Castellated Foam Cushion





Each castellation delivers enhanced support and pressure area care.

Reduces shear and friction forces

The castellations in the *Invacare Propad Profile cushion* reduce shear and friction forces on tissue.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see the [Invacare website](#).

	Size range	Max. user weight	Max. cleaning temperature	Warranty	
Propad Profile Cushion					Manufactured to European Standards: Interior Foam: BS EN 597-1 Cig, BS EN597-2 Match, BS 5852 (Crib 5) Cover: BS 7175 (Crib 5)
	A range of 102 sizes	95 kg (15 stone)	80°	1 year	

Colours

Blue PU

Black Smooth



Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail: ordersuk@invacare.com
www.invacare.co.uk

Invacare Limited
Pencoed Technology Park
Pencoed - Bridgend CF35
5AQ - United Kingdom

Tel. +353 1 8107084
Fax +353 1 8107085
Email: ireland@invacare.com
www.invacare.ie

Invacare Ireland Ltd
Unit 5 - Seatown
Business Campus
Seatown Road
SWORDS - County
Dublin - Ireland



Yes, you can.®

