

AQUATEC®

# Pico Commode



The modular and versatile commode



Yes, you can.®

## Flexible toilet pan holder

You will be surprised how versatile a commode can be when you first see the *Aquatec Pico commode*.

With its modular concept, it can be easily adapted to one's personal demands without the use of any tools. Light weight and completely corrosion free, it sets the new standard for commode chairs.

Choose easily how you would like to handle the toilet pan according to the individual living situation.

Whether you prefer to take out the pan from the rear or rather from the side, the flexible pan holder of the **Pico commode** makes your choice a simple one.

## Support where ever you might need it

Taking off the backrest takes just seconds and the *Aquatec Pico commode* can be used over every toilet for additional support.

## Features and Options



### Versatility

Build your own personal **Pico** without using tools. For limited space e.g. you can use one armrest only.



### High load capacity

High load capacity of 160 kg in every configuration.



## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat height	Seat width	Seat depth	Total width	Total depth	Width between armrests	Weight shower stool	Weight shower chair
<b>Pico Commode</b>	425 - 575 mm	430 mm	420 mm	575 mm	520 mm	430 mm	3,2 / 4,1* kg	4 / 4,8 <sup>1</sup> kg

Max. user weight




**Pico Commode**

160 kg

1. weight with armrests

## Colour

White



Tel: +41 61 487 70 70  
hqeurope@invacare.com  
www.invacare.eu.com

Invacare International GmbH  
Neuhofweg 51  
4147 Aesch BL  
Switzerland



**Yes, you can.®**



© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU - 07/23