

AQUATEC®

Pico



The modular shower chair system



Yes, you can.®

Due to its modular concept and its height adjustability the *Aquatec Pico* can easily be adapted for one's personal needs. The assembly of the light-weight, rustproof shower chair can be realized quick and easy without any tools. Furthermore the product impresses people with its exceptional load capacity.

High loading capacity

With a maximum load capacity of 160 kg the *Aquatec Pico* sets a standard.

Hygiene recess - make your choice

Large hygiene recess either at front or back due to variable position of the backrest.



Features and Options











Aquatec Pico
Shower stool




Aquatec Pico
with armrests

Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat height	Seat width	Seat depth	Total width	Total depth	Width between armrests	Weight shower stool	Weight shower chair
<i>Pico</i>	 425 - 575 mm	 430 mm	 420 mm	 575 mm	 520 mm	 430 mm	 3,2 / 4,1* kg	 4 / 4,8 ¹ kg

Max. user weight	
 160 kg	1. weight with armrests

Colour

White



Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland



© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU - 07/23

Yes, you can.®