

AQUATEC®

Pico



The modular shower chair system



Due to its modular concept and its height adjustability the *Aquatec Pico* can easily be adapted for one's personal needs. The assembly of the light-weight, rustproof shower chair can be realized quick and easy without any tools. Furthermore the product impresses people with its exceptional load capacity.

### High loading capacity

With a maximum load capacity of 160 kg the *Aquatec Pico* sets a standard.

### Hygiene recess - make your choice

Large hygiene recess either at front or back due to variable position of the backrest.



## Features and Options



*Aquatec Pico*

Shower stool



*Aquatec Pico*

with armrests



*Aquatec Pico*

with backrest

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

|             | Seat height  | Seat width | Seat depth | Total width | Total depth | Width between armrests | Weight shower stool | Weight shower chair     |
|-------------|--------------|------------|------------|-------------|-------------|------------------------|---------------------|-------------------------|
|             |              |            |            |             |             |                        |                     |                         |
| <b>Pico</b> | 425 - 575 mm | 430 mm     | 420 mm     | 575 mm      | 520 mm      | 430 mm                 | 3,2 / 4,1* kg       | 4 / 4,8 <sup>1</sup> kg |




Max. user weight

1. weight with armrests

**Pico** 160 kg

### Colour

White



QR Code

Invacare Ltd  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1 656 776 222  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS  
County Dublin - Ireland  
Tel: +353 1 8107084  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

